

Product Data Sheet

Surgical masks




Surgical mask Multiguard antifog
REF 2865 (Vnr. 940001)

Page 1 (2)

Version: 09.08.2013, rev 4

Printed: 5.5.2014

REF	Colour/size	Fixation	
2865	Green, width 17,5 cm	Tie-band	
<p>Classified as type II in the EN14683:2005 European standard Surgical masks – requirements and test methods, Annex A (informative).</p> <p>Bacterial Filtration Efficiency (BFE) >99 %</p> <p>[NATURAL LATEX FREE]</p>			
<p>Raw materials: Outer coverstock Filter Inner coverstock Nose wire Ties Front antifog Back antifog</p>		<p>spunbond polypropylene, blue meltblown polypropylene, white hydrophilic polypropylene, white pe-coated galvanized wire spunbond polypropylene, white laminate white spunbonded polypropylene/transparent polyethylene laminate light blue PPU foam/spunbonded polypropylene</p>	
<p>REACH regulations:</p>		<p>The product covered by this data sheet does not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).</p>	
<p>Features:</p>		<p>An allround antifog surgical mask with high filtration capacity and low breathing resistance. Has foam-edge for antifog functionality and can be used also when wearing glasses. Recommended to be used during dry surgery, with no risk of splashes.</p>	
<p>Pack count:</p>		<p>50 pcs/dispenser, 6 dispensers of 50 pcs (= 300 pcs)/carton</p>	
<p>Dispenser:</p>		<p>Cardboard 185 x 110 x 130 mm</p>	
<p>Transport carton:</p>		<p>Cardboard 280 x 197 x 345 mm</p>	
<p>Sales unit:</p>		<p>Transport carton (6 dispensers of 50 pcs)</p>	
<p>Compliance with legal requirements:</p>		<p>Medical Devices Directive 93/42/EEC, Class I</p>	
<p>Compliance with product standards:</p>		<p>EN 14683:2005, Type II</p>	
<p>Quality standards:</p>		<p>The product is manufactured in compliance with ISO 13485</p>	
<p>Label information:</p>		<p>in accordance with EN ISO 15223-1, EN 980 and EN 1041</p>	
<p>Storage recommendation:</p>		<p>Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sun light and other intensive light sources.</p>	
<p>Disposal recommendation</p>		<p>Follow local instructions</p>	
<p>Risk if reused:</p>		<p>The product is single use. If reused performance of product may deteriorate, cross contamination may occur</p>	

Product Data Sheet

Surgical masks





Surgical mask Multiguard antifog
REF 2865 (Vnr. 940001)

Page 2 (2)

Version: 09.08.2013, rev 4

Printed: 5.5.2014

Product	Dispenser	Transport carton
REF 2865	6438129312973 2865 dispenser  6 438129 312973	6438129512977 2865 carton  6 438129 512977