

PRODUCT SPECIFICATION CATHETER MEDICAL DEVICES

PRODUCT: Curan Advantage Male Catheter

PRODUCT DETAILS:

Code	Size	OD mm	Tip Shape	Colour code
WHM08	08	2.7	Straight	Blue
WHM10	10	3.3	Straight	Black
WHM12	12	4.0	Straight	White
WHM14	14	4.7	Straight	Green
WHM16	16	5.3	Straight	Orange
WHM18	18	6.0	Straight	Red
WHM12C	Ch12	4.0	Tiemann	White
WHM14C	Ch14	4.7	Tiemann	Green

APPEARANCE:

The Advantage Male Catheter is a flexible, DEHP free, PVC, catheter which is supplied sterile. Sterilisation is achieved using irradiation. Each catheter is coated with hydrophilic coating which can be activated using the sachet of 18ml of sterile water included in the pack.

Each device is supplied individually packed in a blister pouch and each pouch contains the catheter with flexible connector and non-touch Blue Grip® for fully hygienic insertion. The sterile water sachet can be activated when the pouch is closed, enabling use of the catheter immediately upon opening the pack.

SPECIFICATION:

Catheter length	42cm
Number of Eyelets	2
Tip	Closed/Straight – all sizes Tiemann – Ch sizes 12 and 14

STORAGE:

Product is supplied sterile so must be stored as received in blister packs and boxes and not opened until use. Do not use if pack is torn or damaged.

USE:

Single use product and labelled as such. Consult instructions for use before using.

SHELF LIFE

2,5 years.

PACKAGING:

Pack Type	Contents
Single Blister Pack	1 Sterile Advantage Hydrophilic Nelaton Catheter - Male
	1 non-touch Blue Grip®
	1 Flexible Connector
	1 Sterile 18ml Water Sachet
Sales Box - Ivory Board	30 Sterile Advantage Hydrophilic Nelaton Catheters – Male
	1 Instructions for Use Leaflet
Outer Transport Carton – Double Corrugated Cardboard	10 Sales Boxes of 30 Sterile Advantage Hydrophilic Nelaton Catheters
	Total Catheters Quantity per Outer -300

PRODUCT PHOTOGRAPH:

