

PRODUCT SPECIFICATION CATHETER MEDICAL DEVICES

PRODUCT: Curan Advantage Female Catheter

PRODUCT DETAILS:

Code	Size	OD mm	Tip Shape	Colour code
WHF08	08	2.7	Straight	Blue
WHF10	10	3.3	Straight	Black
WHF12	12	4.0	Straight	White
WHF14	14	4.7	Straight	Green
WHF16	16	5.3	Straight	Red

APPEARANCE:

The Advantage Female Catheter is a flexible, DEHP free, PVC, catheter which is supplied sterile. Sterilisation is achieved using irradiation. Each catheter is coated with hydrophilic coating which can be activated using the sachet of 8ml of sterile water included in the pack.

Each device is supplied individually packed in a blister pouch and each pouch contains the catheter with flexible connector. The sterile water sachet can be activated when the pouch is closed, enabling use of the catheter immediately upon opening the pack.

SPECIFICATION:

Catheter length 14cm
 Number of Eyelets 2
 Tip Closed/Straight – all sizes

STORAGE:

Product is supplied sterile so must be stored as received in blister packs and boxes and not opened until use. Do not use if pack is torn or damaged.

USE:

Single use product and labelled as such. Consult instructions for use before using.

SHELF LIFE

2,5 years.

PACKAGING:

Pack Type	Contents
Single Blister Pack	1 Sterile Advantage Hydrophilic Nelaton Catheter - Female
	1 Flexible Connector
	1 Sterile 8ml Water Sachet
Sales Box - Ivory Board	30 Sterile Advantage Hydrophilic Nelaton Catheters – Female 1 Instructions for Use Leaflet
Outer Transport Carton – Double Corrugated Cardboard	10 Sales Boxes of 30 Sterile Advantage Hydrophilic Nelaton Catheters Total Catheters Quantity per Outer -300

PRODUCT PHOTOGRAPH:

