

KD-JECT II SYRINGE

EO sterilized and transparent 2-component syringe. The blisterpacked syringe are both free from latex and PVC. KD Ject syringe has smooth and even gliding characteristics.

Properties

- Highly visible scale
- Large wings for easy and firm grip
- Piston stop that prevents the plunger being pushed out
- Color coded dispenser packaging which can be opened in two different ways



802205



805206



805206C



810200



820209

REF	Size	Dispenser box	Transport carton
802205	2 ml Luer Slip centric	100	2500
805206	5 ml Luer Slip eccentric	100	2000
805206C	5 ml Luer Slip centric	100	2000
810200	10 ml Luer Slip eccentric	100	1200
820209	20 ml Luer Slip eccentric	100	800

Raw material

Barrel: Ethylene Polypropylene Copolymer (PP)
Plunger: Polyethylene (PE)

Packaging material

Medical grade paper
Transparent blister packaging PE/PA

Sterilization method:

Ethylene oxide (EO)

Shelf life: 5 years

Medical device directive:

93/94/EEC Class II

Standards:

EN1707, ISO 7886, EN 20594