

Qufora[®] Accessories

Straps | Bed Hanger | Floor stand | Single Packed

Easy-to-use. Durable.

The **Qufora[®] accessories** are practical and durable devices designed for easy use.

PRODUCT BENEFITS

- ▶ The **Qufora[®] leg straps** are designed for the Qufora[®] leg bags. The straps (with two buttons and silicone anti-slip) are adjustable for individual needs. They provide a discrete solution for wearing a urine bag under clothing.
- ▶ The **Qufora[®] hangers** are easy to hang on bed frames.
- ▶ The **Qufora[®] floor stand** comes with four small hooks, making it easy to hang a bed bag on it.
- ▶ The **Qufora[®] accessories** are **single packed** with the REF number, lot number and expiration year/month printed on the protective poly bag.

Straps



Hanger



Floor stand



Qufora® Accessories

Straps | Bed Hanger | Floorstand | Single Packed

Easy-to-use. Durable.

- ▶ **All materials are 100% latex- and DEHP-free.**

QUFORA® ACCESSORIES FEATURES

- The Qufora® **leg straps** come in a small package of two straps (60 cm long). The Qufora® leg straps have a strong holding force. There are silicone strips on the back of the wide bands that ensure correct and reliable positioning on the leg. The leg straps are washable at 40°C.
- With the help of the easy-to-grip tab, the Qufora® **leg straps** can easily be opened even with limited finger dexterity. A strong feature of the Velcro® technology is that it prevents the strap from sticking to other materials, such as stockings.
- The Qufora® **hangers** are durable and fit all beds. The hanger can be cleaned/disinfected in a desk disinfectant for 3 minutes at 85°C. Regular wash in soap and water; combine with wiping with alcohol.
- The Qufora® **floor stand** can be folded flat when not in use.

MATERIALS

- | | |
|---------------|---------------------------------|
| • Straps | Elastic strap: Nylon + Elastane |
| Buttons | ABS |
| Anti-slip | Silicone |
| Velcro | Nylon |
| Wash label | PE |
| Poly bag | PP |
| • Bed hanger | PP/POM |
| Poly bag | PP |
| • Floor stand | PP |

PACKAGING

- | | |
|-----------------------|---|
| • Retail box (Straps) | Cardboard |
| • Retail box (Hanger) | Cardboard |
| • Shipper box | Corrugated cardboard |
| • IFU | None |
| • Shelf-life | Five (5) years |
| • Storage | Dry in temperatures between 0°C and 35°C; protected from direct sunlight; humidity level max 80%; not to be exposed to heavy pressure |

QUFORA® LEG STRAPS

REF	Description	Width / length	Quantity retail box	Quantity shipper box	GTIN-14 ¹
21100301	Velcro® leg straps 1 pair with silicone anti-slip and 2 buttons/strap	3 cm x 60 cm	10 pairs	100 pairs	15710187001715

QUFORA® BED HANGER

REF	Description	Single packed	Quantity retail box	Quantity shipper box	GTIN-14 ¹
21100101	Bed hanger - fits most beds	Yes	5	50	15710187001777

QUFORA® FLOOR STAND

REF	Description	Single packed	Quantity per pack	Quantity shipper box	GTIN-14 ¹
99910003	Floor stand	Yes	1	20	15710187013589

¹GTIN barcode number for retail box

PRODUCT CLASSIFICATION

- Medical Devices conforming to MDR 2017/745 in Class I

APPLIED STANDARDS

- ISO 13485 Quality Management System
- GMDN Straps: 38861; Bed hanger: 41052; Floor stand 41052

Importer

MBH-International A/S
Gydevang 28-30
DK-3450 Allerød
Denmark

Tel: +45 88 30 05 92
info@qufora.com
www.qufora.com



Qufora® is a registered trademark of MBH-International A/S. All rights reserved MBH-International A/S, 3450 Allerød, Denmark.