


GentleCath™ Glide Urinary Catheter with Water Sachet



Intended Use:

The GentleCath™ Glide intermittent catheter is a hydrophilic urinary catheter intended for routine drainage of the bladder. The catheter is designed for transient use and is removed immediately upon emptying the bladder. The product is intended for single use only.

List of Product Range:

CH Size	Outer Diameter (mm)	Male Straight Tip Catheter	Male Tiemann tip catheter	Female Straight Tip Catheter
		Length with funnel: 400 mm Length without funnel: 360 mm	Length with funnel: 400 mm Length without funnel: 360 mm	Length with funnel: 180 mm Length without funnel: 140 mm
		REF	REF	REF
08	2.7	421564	421907	421570
10	3.3	421565	421908	421571
12	4.0	421566	421909	421572
14	4.7	421567	421910	421573
16	5.3	421568	421911	421574
18	6.0	421569	421912	

Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Sterilization:

Sterile	Yes
Sterilized by	Ethylene Oxide
Water Sterilized by	E-beam (beta irradiation)
Shelf Life	18 Months

Certification:

CE-Mark	Yes (CE 0086)
Product Classification	Class I sterile (EU)
Quality System Certification	ISO 13485
Global Medical Device Nomenclature (GMDN)	45603

Material Specification:

Contains Latex	No
Contains DEHP	No
Catheter	Polyolefin-based Elastomer (POBE)
Handling Sleeve	Low Density Polyethylene
Adhesive Tab	Adhesive Film
Sachet Contents	Sterile Purified Water

List of Product Range:								
REF	Product	CH size	Shelf Box		Shipping Carton			MOQ (pcs)
			Qty (pcs)	Dimensions (mm)	Qty (pcs)	Dimensions (mm)	Volume (m3)	
421570	Female Straight Tip	08	30	145 x 93 x 243	240	395 x 309 x 249	0.030	240
421571		10	30	145 x 93 x 243	240	395 x 309 x 249	0.030	240
421572		12	30	145 x 93 x 243	240	395 x 309 x 249	0.030	240
421573		14	30	145 x 93 x 243	240	395 x 309 x 249	0.030	240
421574		16	30	145 x 93 x 243	240	395 x 309 x 249	0.030	240
421564	Male Straight Tip	08	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240
421565		10	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240
421566		12	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240
421567		14	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240
421568		16	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240
421569		18	30	145 x 93 x 446	240	391 x 304 x 455	0.054	240

Manufacturer: **ConvaTec Limited**



Do not use if package has been opened or damaged

For further information on these products, please contact your local ConvaTec sales representative