


100% Silicone Foley Catheter

Standard, 2-Way, Single Packed

Intended use:

The 100% Silicone Foley Catheters are indicated for transurethral drainage of urine from the bladder. The adult 2-way catheters can be retained in the human body for up to 90 days. The standard 2-way catheters (sizes 12 to 18) can also be used for replacement of suprapubic catheters for drainage of the bladder (exchange only). The products are intended for single use only.



Product range:

Silicone Foley Catheter, Standard, 2-Way, Single Packed

| REF* | Total Length mm (±10) | Balloon Capacity ml | CH/FG | Tip | For Suprapubic Exchange | Packaging |
|-------------|-----------------------|---------------------|-------|---------|-------------------------|-----------|
| UM41151210S | 400 | 10 | 12 | Nelaton | Yes | 1/30 |
| UM41151410S | 400 | 10 | 14 | Nelaton | Yes | 1/30 |
| UM41151610S | 400 | 10 | 16 | Nelaton | Yes | 1/30 |
| UM41151810S | 400 | 10 | 18 | Nelaton | Yes | 1/30 |
| UM41152010S | 400 | 10 | 20 | Nelaton | - | 1/30 |
| UM41152210S | 400 | 10 | 22 | Nelaton | - | 1/30 |
| UM41152410S | 400 | 10 | 24 | Nelaton | - | 1/30 |
| UM41152610S | 400 | 10 | 26 | Nelaton | - | 1/30 |



* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Catheter shaft and funnel Silicone
X-ray line Silicone
Balloon Silicone
Tip Silicone, Barium sulphate
Hard valve Polyvinyl chloride, Polypropylene, Acrylonitrile-butadiene-styrene, Silicone, Stainless steel
Sleeve Polypropylene

Certifications:

CE-marking? Yes (0123)
Product classification Class IIB (EU)
Quality system certified to ISO 13485
Global Medical Device Nomenclature (GMDN)34917
 For suprapubic exchange (sizes 12 to 18): 34924

Sterilization:

Sterile? Yes
Sterilized by Ethylene oxide
Shelf life 5 years

100% Silicone Foley Catheter

Standard, 2-Way, Single Packed



Packaging specifications:

| REF | Primary package/ Single peelpack | | Secondary package/ Shelf box | | Tertiary package/ Shipper | | | | Pallet |
|-------------|-------------------------------------|-------------------------|---------------------------------|---------------------------|------------------------------|---------------------------|-----------------|--------------------------|----------|
| | Qty, pcs | Dimensions, mm (LxW) | Qty, pcs | Dimensions, mm (LxWxD) | Qty, pcs | Dimensions, mm (LxWxD) | Gross Wt. kg | Volume m ³ | Qty, pcs |
| UM41151210S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41151410S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41151610S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41151810S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41152010S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41152210S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41152410S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |
| UM41152610S | 1 | 500 x 66 | 1 | 510 x 80 x 20 | 30 | 525 x 225 x 260 | 3.0 | 0.031 | 720 |





EAN codes:

| REF | Primary package | Secondary package | Tertiary package |
|-------------|--------------------|--------------------|--------------------|
| UM41151210S | (01)30768455153856 | (01)50768455153850 | (01)70768455153854 |
| UM41151410S | (01)30768455153863 | (01)50768455153867 | (01)70768455153861 |
| UM41151610S | (01)30768455153870 | (01)50768455153874 | (01)70768455153878 |
| UM41151810S | (01)30768455153887 | (01)50768455153881 | (01)70768455153885 |
| UM41152010S | (01)30768455153894 | (01)50768455153898 | (01)70768455153892 |
| UM41152210S | (01)30768455153900 | (01)50768455153904 | (01)70768455153908 |
| UM41152410S | (01)30768455153917 | (01)50768455153911 | (01)70768455153915 |
| UM41152610S | (01)30768455153924 | (01)50768455153928 | (01)70768455153922 |

Packaging material:

| | |
|-------------------|---|
| Primary package | Sleeve - Ethylene-vinyl acetate, Peelpack pouch - Coated paper, Polyethylene terephthalate/Polyethylene |
| Secondary package | Ivory cardboard |
| Tertiary package | Corrugated board |

Storage and handling:

| | |
|---|---|
|  | FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care |
|  | KEEP DRY |
|  | THIS WAY UP - Indicates correct upright position of the transport package |
|  | DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package |


Well Lead Medical Co., Ltd.
 C-4 Jinhu Industrial Estate
 Hualong, Panyu, Guangzhou
 511434 P.R. China

Made in China

EC REP

Shanghai Intl. Holding Corp. GmbH
 Eiffestrasse 80, 20537 Hamburg
 Germany

Distributed by



ConvaTec Limited
 First Avenue, Deeside Industrial Park
 Deeside, Flintshire, CH5 2NU, UK
www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-019441-MM, ver. OCT 2018