



evercare[®] Medi Blanket

REF: 5300, 5301

evercare[®] Medi Blanket is available in two sizes, choose between full length and half-length for thermal comfort and warmth. Single use blanket for patients in healthcare settings, manufactured in controlled environment to support infection control.



Key Features

- ▶ Single use
- ▶ Available in two sizes

REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
5300	Medi Blanket, single pack, non-sterile, White	Large 130x210 cm	1 / 34
5301	Medi Blanket, single pack, non-sterile, White	Small 130x110 cm	1 / 40

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Main material	3 layers: PP - 70% PES and 30% Viscose - PP

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	PE bag, Transparent	500 x 350
Dispenser Pack	PE bag, Transparent	5300: 590 x 360 x 825 5301: 560 x 470 x 730
Transport Carton	Carton, TCF or ECF quality	5300: 580 x 350 x 550 5301: 580 x 350 x 390
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
5300	6438129122299	N/A	6438129522297
5301	6438129122305	N/A	6438129522303

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.