

## DOCTOR'S

EB00.003

### Doctor's bag

The design of this bag is ideal for healthcare professionals that make home visits, because it allows them to transport all the medical instruments they may need. Its large capacity, internal layout and ergonomic features make this briefcase an essential tool for doctors that perform their activity outside the medical center.



#### Description

- Front pocket for documents & clipboard
  - Organizer pocket in the interior
  - Hidden backpack straps
  - Trolley band and External ID badge holder
  - The large outer file-folder pocket padded and will accommodate up to a 12" laptop
  - 2 net pockets for instruments
  - Main compartment with 4 removable pockets with transparent window
  - Ergonomic padded handle
  - Zipper for padlock
  - Isothermal ampoule holder with 38 elastic bands of various sizes
  - Large flat pocket for documents
- [embed]https://www.youtube.com/watch?v=l-pFzKjwYKg[/embed]

#### Additional information

Weight	2.6 Kg
Dimensions	17 x 40 x 30 cm
Color	BLACK
Capacidad (Litros)	20.40
Contenido máximo recomendado (Kg)	5.00