


# Purifeed

## Feeding Tubes

**Intended use:**

A feeding tube is inserted through the nose or mouth via oesophagus into the stomach, in order to administer medicine, liquids and/or enteral feeding by means of a syringe or tubing. The product is intended for single use only.



**List of product range:**

Purifeed feeding tubes	Effective length mm				
		Tube, Outer diameter, mm			
		1.67	2.0	2.67	3.33
		CH05 REF*	CH06 REF*	CH08 REF*	CH10 REF*
	375	3100118x	3100218x	3100318x	3100418x
	475	3100518x	3100618x	3100718x	
<b>X-ray line</b>	975	3100818x	3100918x	3101018x	3101118x

All REF: Shore A hardness : 78

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

**Material specifications:**

**Contains latex?** . . . . . No  
**Contains DEHP?** . . . . . No  
**Tube** . . . . . Polyurethane with X-ray line  
**Connector** . . . . . Polyvinyl chloride

**Certification:**

**CE-marking?** . . . . . Yes (CE 0086)  
**Product classification** . . . . . Class IIB (EU)  
**Quality system certified to** . . . . . ISO 13485  
**Global Medical Device Nomenclature (GMDN)**. . . . . 14221

**Sterilisation:**

**Sterile?** . . . . . Yes  
**Sterilized by** . . . . . Ethylene oxide  
**Shelf life** . . . . . 5 years

# Purifeed

## Feeding Tubes

### Packaging:

REF*	Shelf box		Shipping carton				MOQ, pcs.
	Qty., pcs.	Dimensions, mm	Qty., pcs.	Dimensions, mm	Gross Wt., kg	Volume, m <sup>3</sup>	
3100118x	100	145 x 74 x 576	1000	391 x 304 x 585	11.0	0.070	1,000
3100218x	100	145 x 74 x 576	1000	391 x 304 x 585	11.0	0.070	1,000
3100318x	100	145 x 74 x 576	1000	391 x 304 x 585	11.0	0.070	1,000
3100418x	100	145 x 74 x 576	1000	391 x 304 x 585	11.0	0.070	2,000
3100518x	100	145 x 74 x 734	1000	391 x 304 x 743	11.0	0.088	3,000
3100618x	100	145 x 74 x 734	1000	391 x 304 x 743	11.0	0.088	3,000
3100718x	100	145 x 74 x 734	1000	391 x 304 x 743	11.0	0.088	3,000
3100818x	100	294 x 125 x 360	600	391 x 304 x 743	13.5	0.088	2,400
3100918x	100	294 x 125 x 360	600	391 x 304 x 743	13.5	0.088	1,800
3101018x	100	294 x 125 x 360	600	391 x 304 x 743	13.5	0.088	1,800
3101118x	100	294 x 125 x 360	600	391 x 304 x 743	13.5	0.088	2,400

\* Here noted as x, reference codes contain an additional digit (1-5) representing a language code.  
 Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

### Storage and handling:



*With reservations for misprints and alterations in product specifications. Printed version not controlled.*

Manufacturer: **ConvaTec Limited**  
 First Avenue  
 Address: Deeside Industrial Park  
 Deeside, Flintshire, CH5 2NU, UK

For further information on these products, please contact your local ConvaTec sales representative.

TM indicates a trademark of Unomedical a/s.  
 The ConvaTec logo is a trademark of ConvaTec Inc.

