

4232 Medical Face Mask Extra Protection

Medical Face Mask

Product details

Size : One size
Descriptive feature : Anti-fog, Tie Band, Type IIR, Visor
Colour : Blue
Sterile : Non-sterile
Brand : BARRIER

Images



Delivered items

4232-01

Country of origin: Japan

Shelf life: 5 years

Sterilization method: Non-sterile

Packing information: First packaging layer is a paper board dispenser box. Second layer is a corrugated board transport box.

Is suitable for Tray: Yes

Packing level	Quantity	GS1 code	WxLxH (mm)	Vol (dm3)	Weight gross/net (kg)
Dispenser box	50	7323190043255	200x310x105		
Transport box	200	7323190043248	208x318x436	28,8	3.1 / 2.1
Pallet	8000	7323190043231	800x1200x1894	1818.2	

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, Gamlestadsvägen 3 C, SE-402 52 Göteborg, Sweden. Phone +46 31 722 30 00.
 The Mölnlycke trademarks, names and logotypes are registered globally to one or more of the Mölnlycke Health Care Group of Companies.

©2021 Mölnlycke Health Care AB. All rights reserved.



Material

Animal tissues :	No
Human blood derivatives :	No
Natural rubber latex :	No
Medicinal substances :	No
Phthalates :	No
Polyvinyl chloride :	No

Product Composition Medical Face Masks

Product Component	Composition
Outer facing	Cellulose/Polyester
Inner facing	Cellulose/Polyester
Filter medium	Polypropylene nonwoven
Insertion	Polypropylene nonwoven
Nose clip	Aluminium
Tie bands	Polyester
Visor	Polyester
Anti-fog film	Polypropylene nonwoven coated with polyethylene film

Product Performance Medical Face Masks, Type IIR

Characteristics	Test Method	Internal Test Method	Unit	Requirement	Product Performance
Bacterial filtration efficiency (BFE)	EN 14683 (Annex B)	N/A	%	≥98	≥ 99
Differential pressure	EN 14683 (Annex C)	N/A	Pa/cm ²	<60	26-36
Splash resistance pressure	ISO 22609	N/A	kPa	≥16.0 Less than 3 penetrations out of 32 tested	31 of 32 pass at 16 kPa
Microbial cleanliness	ISO 11737-1	T-303	CFU/g	≤30	0-10

Find out more at www.molnlycke.com

