

evercare[®] Adhesive OP-tape, sterile

REF: 1493-01

evercare[®] Adhesive OP-tape, sterile is a handy component for fixation of drapes, tubings etc. during surgical procedures.



10 x 58 cm

Key Features

- ▶ Fixation of drapes, tubings etc
- ▶ Transparent


REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
1493-01	Adhesive OP-tape, sterile, Transparent	10 x 58 cm	10 / 200

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-07 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN ISO 13795, EN ISO 11135-1, EN ISO 11737-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Main material	Polyurethane film coated with pressure sensitive acrylic adhesive

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	Paper, PE/PET film	350 x 415
Dispenser Pack	Carton, TCF or ECF quality	381 x 279 x 320
Transport Carton	Carton, TCF or ECF quality	574 x 401 x 330
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
1493-01	N/A	6438129318117	6438129518111

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.