

SELEFA® Optimia Nonwoven swabs

REF: 223001, 223002, 223003, 223004

SELEFA® Optimia Nonwoven swabs, 4-ply are produced from soft and highly absorptive nonwoven material. They are primarily used for cleaning and exudate absorption of wounds and abrasions.

Key Features

- ▶ High quality viscose/polyester nonwoven material
- ▶ Low linting and high softness
- ▶ Keeps its shape, even when wet
- ▶ Free of binding agents and optical brighteners
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than gauze devices
- ▶ Available in several sizes
- ▶ For external wound care, treating burns or abrasions, in the first aid area etc.
- ▶ Can be sterilized using EO
- ▶ Optimized packaging: convenient dispensing and better protection from contamination



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
223001	Optimia Nonwoven swabs, 4-ply, non-sterile	5x5 cm	150 / 13500
223002	Optimia Nonwoven swabs, 4-ply, non-sterile	7,5x7,5 cm	150 / 9000
223003	Optimia Nonwoven swabs, 4-ply, non-sterile	10x10 cm	150 / 5400
223004	Optimia Nonwoven swabs, 4-ply, non-sterile	10x20 cm	75 / 900

Sales Unit: Transport Carton