

# SELEFA® PE gloves, Standard, X-Long

REF: 210362M

PE gloves Standard, X-Long, Transparent, embossed waterproof personal protective gloves used primarily in following applications: Salons / spas, electronic assembly, packaging, automotive, janitorial, general pharmaceutical and home use.



## Key Features

- ▶ Transparent, embossed
- ▶ Made of polyethylene

## Recommendations for use



Basic cleaning



Wash

REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
210362M	PE gloves, Standard, X-Long, embossed, M, Transparent	■	530 mm	100 / 2000

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for personal protective equipments.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Compliance with legal requirements:</b>	Personal Protective Equipment Regulation (EU) 2016/425, PPE R category III EU Declaration of Conformity for Personal Protective Equipment: <a href="https://www.evercare-medical.com/about-us">https://www.evercare-medical.com/about-us</a>
<b>Food safety:</b>	In accordance with EU Regulation 1935/2004 and 10/2011 the glove is suitable for handling all foodstuffs.
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
Palm and Finger	LDPE & LLDPE, single wall thickness min. 0.025 mm

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Cardboard, TCF or ECF quality	230 x 135 x 35
<b>Transport Carton</b>	Corrugated board, TCF or ECF quality	370 x 245 x 290
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>210362M</b>	N/A	6438129354638	6438129554632

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.