

# SELEFA® Round plaster, nonwoven

REF: 240102F, 240103F

SELEFA® Round plaster, nonwoven is an injection plaster with lined pad to prevent sticking to the wound.



## Key Features

- ▶ Lined pad to prevent sticking to the wound


REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
240102F	Round plaster, nonwoven, Beige	3,8 x 3,8 cm	200 / 10200
240103F	Round plaster, nonwoven, Beige	2 x 4 cm	200 / 10200

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	ASTM D3330
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>Carrier</b>	Nonwoven, 100 % PES
<b>Pad, surface</b>	LDPE
<b>Pad, absorbent core</b>	Cellulose/PP/PES
<b>Skin adhesive</b>	Synthetic, colofonium-free
<b>Release paper</b>	Silicone

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	240102F: 92 x 165 x 35 240103F: 93 x 165 x 35
<b>Transport Carton</b>	Carton, TCF or ECF quality	580 x 185 x 300
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>240102F</b>	N/A	6438129352191	6438129552195
<b>240103F</b>	N/A	6438129352207	6438129552201

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.