

# Terumo

## Dental Needle



On this page: [Description](#) [Item specifications](#)

### Description

Terumo Dental Needles are designed with the goal to minimize patient trauma.

The dental needles are available in short and long sizes and with imperial and metric thread.

The needles are suitable for both conduction anaesthesia and infiltration anaesthesia.

---

### Item specifications

Reference <sup>1</sup>	Cannula			Hub	
	Gauge (G)	Effective Length (mm)	Length at cartridge (mm)	Type	Colour
<b>DN*M2508F</b>	25	8	10.7	Metric	Orange
<b>DN*M2535F</b>	25	35	10.7	Metric	Orange
<b>DN*M2541</b>	25	41	10.7	Inch/Imperial	Orange
<b>DN*M2716F</b>	27	16	10.7	Metric	Grey
<b>DN*M2721</b>	27	21	10.7	Inch/Imperial	Grey
<b>DN*M2721B</b>	27	21	13.2	Metric	Grey
<b>DN*M2721F</b>	27	21	10.7	Metric	Grey
<b>DN*M2722</b>	27	22	10.7	Inch/Imperial	Grey
<b>DN*M2730</b>	27	30	10.7	Inch/Imperial	Grey
<b>DN*M2735</b>	27	35	10.7	Inch/Imperial	Grey
<b>DN*M2735B</b>	27	35	13.2	Metric	Grey
<b>DN*M2741</b>	27	41	10.7	Inch/Imperial	Grey
<b>DN*M3013B</b>	30	13	13.2	Metric	Yellow
<b>DN*M3013F</b>	30	13	10.7	Metric	Yellow
<b>DN*M3016</b>	30	16	10.7	Inch/Imperial	Yellow
<b>DN*M3016B</b>	30	16	13.2	Metric	Yellow
<b>DN*M3016F</b>	30	16	10.7	Metric	Yellow
<b>DN*M3019</b>	30	19	10.7	Inch/Imperial	Yellow
<b>DN*M3021</b>	30	21	10.7	Inch/Imperial	Yellow
<b>DN*M3021B</b>	30	21	13.2	Metric	Yellow
<b>DN*M3021F</b>	30	21	10.7	Metric	Yellow
<b>DN*M3022</b>	30	22	10.7	Inch/Imperial	Yellow

 Enlarge table

Unit box: 100 pieces

Shipping carton: 2000 pieces