

VICRYL*
(Polyglactin 910) Braided and Monofilament

SUTURE

Description:	Braided	For easy handling and secure knot tying.
	Coated	For smooth passage through tissue and easy knot tie down.
	Synthetic	For minimal tissue reaction.
	Absorbable	By hydrolysis. Essentially complete in 56-70 days.

Tensile strength:	Post implantation	Approximate % original strength remaining
	14 days:	75%
	21 days:	50% (6-0 and larger) 40% (7-0 and smaller)
	28 days:	25% (6-0 and larger)

Color: Violet or undyed.









Range: 10-0 to 2 (USP). Supplied as needled sutures and ligatures. Sizes 10-0 and 9-0 are monofilament in structure.



Indications: VICRYL[®] sutures are intended for use in soft tissue approximation and/or ligation, including use in ophthalmic surgery, peripheral nerve adaptation and microsurgery for vessels less than 2 mm in diameter. Typical areas of use include:

- Fascia closure - due to the unique 28-day claim
- Subcutaneous fat
- Joint capsule
- Uterus

Full details in the Instruction for Use included in every package.



▼ 3/8 Circle Reverse Cutting PRIME [®] Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 P-1 11 mm	6-0	○	45	MULTIPASS [®]	36	MPV489H
	0-0	○	45		30	V489H
	5-0	○	45	MULTIPASS [®]	36	MPV490H
	5-0	○	45		36	V490H
 P-3 13 mm	6-0	○	45	MULTIPASS [®]	36	MPV492H
	6-0	○	45		36	V492H
	5-0	○	45	MULTIPASS [®]	36	MPV493H
	5-0	○	45		36	V493H
	4-0	○	45	MULTIPASS [®]	36	MPV494H
	4-0	○	45		36	V494H
 FS-3 16 mm	6-0	○	45	MULTIPASS [®]	36	MPV500H
	4-0	○	45		24	W9506T
 PS-2 19 mm	5-0	○	45	MULTIPASS [®]	36	MPV495H
	4-0	○	45	MULTIPASS [®]	36	MPV496H
	4-0	○	40		24	W9515T
	3-0	○	45	MULTIPASS [®]	36	MPV497H
	3-0	○	45		36	V497H
 PS-2S 10 mm	4-0	○	45	Slimmed	36	V498H
 PS-1 24 mm	3-0	○	45	MULTIPASS [®]	36	MPV893H
	3-0	○	40		30	V893H
	2-0	○	45	MULTIPASS [®]	36	MPV894H
	2-0	○	45		36	V894H
 PS 26 mm	3-0	○	70	MULTIPASS [®]	36	MPV9526H
 PSL 30 mm	2-0	●	70		36	V886H

▼ 3/8 Circle Reverse Cutting Needles						
Needle	Gauge	Color	Length (cm)	Other information	Number per box	Code
 C-3 13 mm	6-0	●	45		36	V384H
	5-0	●	45		36	V385H
	4-0	●	45		36	V386H
 FS-3 10 mm	6-0	●	45		36	V387H
	5-0	●	45		36	V388H