

# HeartStart

## Four-Year Battery

M5070A

[Find similar products >](#)

This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.

Contact & support



**Contact us >**

**Request contact**

## Specifications

### Accessories Battery

Technology Lithium

Life Time 4 years

### Product details

Use with Philips Healthcare Equipment M5066A, M5067A, M5068A, 861304

Use with Non-Philips Healthcare Equipment No

[See all specifications >](#)