

22





Sterowash

The 2-in-1 eye and wound wash could save someone's sight!

In an emergency a twist off top makes Sterowash ready for use in a second!

It's also perfect to cleanse wounds to reduce the risk of infection, as its sterile.

Once fixed on the wall where it'll be seen every day, it's there permanently in case of emergency.

With a slimline neck for improved directional accuracy, and a twist off top ready in the blink of an eye, these pods are the only choice in an emergency.

Sterowash saline 20ml pods treat both eyes and wounds and are available on **Drug Tariff**, e.g. as prescribed by your doctor. Class IIA Medical device.

The Ultra Durable Premier Box

Ideal for warehouses, factories, garages etc, as it protects the eye and wound wash, and keeps it clean. These stations are specifically designed for your workplace, you decide which kit you prefer.

Description	ltem	Ref
Sterowash Eyecare Kit 2 (2 x 500ml)	1	9114
Sterowash Eyecare Kit 3 (3 x 500ml)	1	9116
Sterowash 500ml bottles - Shelf life 4 years 112 cases per pallet	1x10	1405



Eye & Woundwash Kit, with wall bracket and 2 x 500ml bottles & 5 x 20ml pods.



Sterowash Eyewash Kit 2, with (2 x 500ml) bottles of eyewash with eyepads provided to cover the eye while waiting for medical attention.