



Stetoskop vuxen

Artikel: Adscope 601 Convertible

Pris enligt avtal: 779,00 st.

Livstid garanti

Beställes via ONEMED

“Combining AFD technology with a convertible design makes this the most versatile of our cardiology scopes”

Producent: American Diagnostic Corporation



Dosmanometer för blodtrycksmätare

Artikel: 808N, Dosmanometer, Gauge F/760

Pris enligt avtal: 102,00 st.

Garanti kalibrering 20 år

Beställes via ONEMED

“Gauge, 300mmHg, Individually tested to conform to ANSI/AAMI/ISO 81060-1 standards, All gauges individually boxed”

Producent: American Diagnostic Corporation



Bollpump med ventil för blodtrycksmätare

Artikel: 875N, bollpump med ventil

Pris enligt avtal: 37,00 st.

1 års garanti

Beställes via ONEMED

“Bulbs for precise control during inflation and deflation insist on ADC’s bulbs and valves”

Producent: American Diagnostic Corporation



Blodtrycksmanschett lårben

Artikel: 845-13TBK-2, blodtrycksmanschett

Pris enligt avtal: 205,00 st.

3 års garanti

Beställes via ONEMED

“Adcuff and Bladders, 2 Tube, Durable nylon construction, 2-tube inflation bladder with 25cm tubing length”

Producent: American Diagnostic Corporation