



optilube.™

with 2% lidocaine
for urethral catheterisation

active CHG free

4g less plastic per syringe ✓

1.7
TONNES
SAVED!

That's the equivalent
of an entire hippopotamus
weight in plastic,
saved every single year.



OptiLube Active CHG Free is a single-use, sterile, pre-filled syringe, containing local anaesthetic (Lidocaine 2%). It now comes packed in a soft pouch, with extra peel area, to reduce plastic and make it even easier to open.

That means we can save almost 2 tonnes of plastic every year by making the switch to the new **OptiLube Active CHG Free**.

The most important step that clinicians can take to help reduce the risk of CAUTI is using an effective aseptic technique during the catheterisation procedure.

OptiLube Active CHG Free, as with each catheterisation jelly in the **OptiLube** range, reduces the risk of pain, trauma and infection during catheterisation by coating the full length of the urethra in lubricant.

OptiLube Active CHG Free also provides added reassurance to clinicians that patients with sensitivities to Chlorhexidine will not experience allergic reactions due to this active ingredient.



CODE	DESCRIPTION	INNER PACK	OUTER PACK
1168S	OptiLube Active CHG Free 11ml pre-filled syringe with soft packaging	Pack of 10	Case of 200

As healthcare lubrication specialists, our high quality yet cost-effective range of catheterisation gels offers clinicians a choice. This means the most appropriate format of gel can be picked for each individual patient.

The range includes...

optilube.™

optilube.™
active

optilube.™
active CHG free

All of our catheterisation gels are:

- Available in 6ml (female) and 11ml (male) volumes
- Provided in a sterile, single-use syringe format
- Developed to achieve the perfect viscosity for lubricating the full length of the urethra, helping to reduce the risk of pain, trauma and infection.