

Product Data Sheet

SOFTCARE® Face Shield

REF MI98602

REF	Colour/size	Fixation	
MI98602	Clear / width 33 cm x height 22 cm	Elastic band	
Raw materials:			
<ul style="list-style-type: none"> - Visor - Forehead foam pad - Elastic band 			<ul style="list-style-type: none"> - Polyethylene terephthalate (PET) film, antifog, clear - Polyurethane foam, white - Synthetic rubber, white
Brand:			Softcare®
Features:			A face shield with clear plastic visor. Comfortable forehead foam pad. Elastic band attachment.
Pack count:			80 pcs/carton, each product individually packed in PE bag
Sales unit:			Carton
Storage recommendation:			Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sun light and other intensive light sources.
Disposal recommendation:			Follow local guidelines
Risk if reused:			The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Product	
REF MI98602	