

Gauze swabs, 5-pack, sterile

REF: 2222705, 2222725, 2222745

SELEFA® Gauze swabs, 5-pack, sterile can be used for non-invasive medical procedures, such as wound padding, absorbing blood and exudates, skin and wound cleansing or as secondary dressing.

Key Features:

- ▶ 100% cotton gauze (17 threads per cm²)
- ▶ Peroxide bleached (no chlorine)
- ▶ 8-ply
- ▶ Folded in cut edges to prevent loose threads
- ▶ Sterilized by ethylene oxide gas (EO)
- ▶ 5 swabs per pouch



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2222705	Gauze swabs, sterile	5 x 5 cm, 8-ply	125/2500
2222725	Gauze swabs, sterile	7,5 x 7,5 cm, 8-ply	125/2500
2222745	Gauze swabs, sterile	10 x 10 cm, 8-ply	125/2250

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices
EO sterilized according to valid version of ISO 11135.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5; EN ISO 10993-7; EN ISO 10993-10

Compliance with legal requirements: Medical Devices Directive 93/42/EEC, Class Is

Compliance with product standards: EN 14079

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 11607-1 and EN 1041

Material

Product Part	Raw Material & Characteristic
Raw material	100 % cotton gauze, white, 17 threads/cm ²

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical paper / Plastic film PA/PE	2222705: 75 x 90
		2222725: 100 x 114
		2222745: 133 x 135
Dispenser Pack	Carton, ECF or TCF quality	2222705: 90 x 78 145
		2222725: 115 x 95 x 145
		2222745: 135 x 130 x 145

Transport Carton	Carton, ECF or TCF quality	2222705: 400 x 190 x 300
		2222725: 600 x 200 x 300
		2222745: 400 x 300 x 420
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2222705	 6 438129 151015	 6 438129 351019	 6 438129 551013
2222725	 6 438129 151022	 6 438129 351026	 6 438129 551020
2222745	 6 438129 151039	 6 438129 351033	 6 438129 551037



Manufactured in China for
OneMed Group Oy

Address: Metsäläntie 20
Postal Address: P.O. Box 10
FI-00321 Helsinki, Finland
Phone +358 20 786 6800
www.onemed.com
Domicile Helsinki
Business ID 2039640-1