


evercare[®] Fluid collection pouch with sieve, 62 x 50 cm

REF: 1873

evercare[®] Fluid collection pouch with sieve, sterile is a transparent irrigation pouch which contains an adhesive strip, sieve and exit port. Designed to collect irrigation fluids and tissue specimens during surgical procedures.



Key Features


- ▶ Transparent
- ▶ Adhesive bag

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1873	Fluid collection pouch with sieve, sterile, Transparent	62 x 50 cm	11 / 44

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-01
EN ISO 10993-7

Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, MDR class Is

Compliance with product standards:

EN ISO 13795, EN ISO 11135:2014/A1:2019

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Main material	100% PEVA, transparent
Sieve	PE/PET nonwoven, white
Valve for fluid collection pouch	PE, white
Wire for fluid collection pouch	Metal, coated with PE

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	250 x 390
Dispenser Pack	TCF or ECF quality	285 x 242 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1873	06438129111750	06438129311754	06438129511758

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.