


evercare[®] Banded bag

REF: 1850, 1855-01

evercare[®] Banded bag, sterile are circular, disposable bags for covering a variety of equipment for use in the operation room and other hospital departments where sterile equipment draping is needed.

Key Features

- ▶ Transparent
- ▶ Circular disposable bag
- ▶ Complies to EN 13795-1 High performance critical product area requirements




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1850	Banded bag, sterile, Transparent	Ø 90 cm	18 / 72
1855-01	Banded bag, sterile, Transparent	Ø 140 cm	10 / 40

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-01
EN ISO 10993-7

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards: EN13795-1

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Bag	PE-film, transparent
Elastic band	Latex free, blue

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, Transparent	150 x 290
Dispenser Pack	TCF or ECF quality	265 x 143 x 234
Transport Carton	TCF or ECF quality	545 x 310 x 246
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1850	06438129111712	06438129311716	06438129511710
1855-01	06438129119107	06438129319101	06438129519105

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.