

SELEFA® Nonwoven round sponges, 5-pack, sterile

REF: 22325030, 22325040, 22325050

SELEFA® Nonwoven round sponges, sterile have a wide area of use: pre-operative skin cleansing, wound and abrasion cleansing, wiping, absorption, protection etc. They are also intended to stop bleeding or to apply topical medications. Their high softness insure lenient contact with wound and skin as the material is less adherent and less linting than gauze devices.

Key Features

- ▶ High quality viscose/polyester nonwoven material
- ▶ High softness and low linting
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than woven gauze devices
- ▶ Sterilized using EO
- ▶ Easy-to-use unit pouch
- ▶ 5-pack




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
22325030	Nonwoven round sponges, 5-pack, sterile	Ø 30 mm	75 / 1500
22325040	Nonwoven round sponges, 5-pack, sterile	Ø 40 mm	75 / 1350
22325050	Nonwoven round sponges, 5-pack, sterile	Ø 50 mm	50 / 900

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards: EN ISO 11607-1, EN ISO 11135

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 70% viscose and 30% polyester, white
Elastic ring	Silicone-based elastic

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	22325030: 98 x 113 22325040: 122 x 138 22325050: 122 x 138
Dispenser Pack	Carton, TCF or ECF quality	22325030: 115 x 95 x 140 22325040: 140 x 120 x 150 22325050: 140 x 120 x 150
Transport Carton	Carton, TCF or ECF quality	22325030: 495 x 240 x 290 22325040: 435 x 370 x 310 22325050: 435 x 370 x 310
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
22325030	06438129151183	06438129351187	06438129551181
22325040	06438129151190	06438129351194	06438129551198
22325050	06438129151206	06438129351200	06438129551204

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.