



SELEFA® Optimia Flat gauze sponges

REF: 222052

SELEFA® Optimia Flat Gauze sponges are highly absorbent sponges for a wide area of use such as: skin cleansing, wound cleansing, wiping, fluid and moisture absorption, protection, disinfection etc. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze. The compact folding keeps the product shape stable. Gauze sponges in these small sizes are suitable for various applications as cleaning of small wounds and abrasions, minor prepping, debris removal etc.

Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ Folded edges prevent loose threads
- ▶ Can be sterilized using EO
- ▶ 17 threads per cm²




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
222052	Optimia Flat gauze sponges, 100-pack, non-sterile, White	30x30 mm	1000 / 10000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 14079
Other tests performed:	0
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	PE bag, Transparent	170 x 210
Dispenser Pack	PE bag, Transparent	270 x 400 x 0
Transport Carton	Carton, TCF or ECF quality	3395 x 290 x 360
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
222052	06438129154245	06438129354249	06438129554243

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.