



SELEFA® Triangular nonwoven sponges, sterile

REF: 223241

SELEFA® Triangular nonwoven sponges, sterile have a wide area of use: pre-operative skin cleansing, wound and abrasion cleansing, wiping, absorption, protection etc. They are also intended to stop bleeding or to apply topical medications. Their high softness insure lenient contact with wound and skin as the material is less adherent and less linting than gauze devices.

Key Features

- ▶ High quality viscose/polyester nonwoven material
- ▶ High softness and low linting
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than woven gauze devices
- ▶ Sterilized using EO
- ▶ Easy-to-use unit pouch
- ▶ 5-pack




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
223241	Triangular nonwoven sponges, 5-pack, sterile, White	30x30x30 mm	250 / 5000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN ISO 11607-1, EN ISO 11135, EN 29073, EN 1644-1,-2
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 70% viscose and 30% polyester, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	
Dispenser Pack	Carton, TCF or ECF quality	237 x 100 x 116
Transport Carton	Carton, TCF or ECF quality	595 x 245 x 415
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
223241	06438129154160	06438129354164	06438129554168

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.