



## evercare® XP Special OP-Gowns, sterile

REF: 2534, 2535, 2537, 2539, 2540, 2544, 2545, 2547, 2549, 2550

evercare® XP Special OP-Gowns are suitable for procedures with medium to high quantities of fluids. Made in a soft SMS material with antistatic and alcohol repellent treatment for enhanced user protection. evercare® XP Special is reinforced with a breathable material in the critical zone, in the front and in the sleeves, to maximize protection against contamination of the surgical wound and to protect the user at wet procedures. All evercare® OP-Gowns provide an effective barrier to body fluids, skin particles and microbes and are designed for optimal comfort and ease of movement. The adjustable Velcro in the neck enables optimal fit while securing adjustments. The gown is sterilized with EO-gas and packed in easy-to-open unit pouches. State-of-the-art packaging system with handy functional dispenser and labeling. Unit packs are supplied with barcoded removable labels for improved traceability.

### Key Features

- ▶ Good protection and Xtra comfort, antistatic and alcohol repellent treated
- ▶ Ultrasonic heat sealed seams for sleeves and shoulders to increase comfort by reducing excessive materials causing discomfort
- ▶ Gown with overlapping back and adjustable Velcro in the neck enables optimal fit for individual adjustments
- ▶ Clear stamp with size and performance type
- ▶ Belt card with clear indications of sterile / non-sterile for secure donning
- ▶ High quality cuffs for improved comfort and tight fit round wrists
- ▶ Book-folding enabling unpacking and donning without risk for contamination
- ▶ SMS material with low linting and very high abrasion resistance




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2534	XP Special OP-Gown, sterile, Blue	M-short, Length 120 cm	21 / 42
2535	XP Special OP-Gown, sterile, Blue	M, Length 131 cm	20 / 40
2537	XP Special OP-Gown, sterile, Blue	L, Length 142 cm	20 / 40
2539	XP Special OP-Gown, sterile, Blue	XL, Length 152 cm	18 / 36
2540	XP Special OP-Gown, sterile, Blue	XXL, Length 160 cm	18 / 36
2544	XP Special OP-Gown with hand towels, sterile, Blue	M-short, Length 120 cm	21 / 42
2545	XP Special OP-Gown with hand towels, sterile, Blue	M, Length 131 cm	20 / 40
2547	XP Special OP-Gown with hand towels, sterile, Blue	L, Length 142 cm	20 / 40
2549	XP Special OP-Gown with hand towels, sterile, Blue	XL, Length 152 cm	18 / 36
2550	XP Special OP-Gown with hand towels, sterile, Blue	XXL, Length 160 cm	18 / 36

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class Is
Compliance with product standards:	EN ISO 11135, EN 13795-1 High performance
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

## Material

Product Part	Raw Material & Characteristics
Body, sleeve, belt	PP SMS nonwoven, blue
Reinforcement sleeve and front	PP nonwoven / breathable PE-film, white
Cuff	PES, rib knitted fabric
Neck ribbon	PP Spunbond nonwoven, white
Velcro band	PA woven ribbon with hook or loop
Belt card	Cellulose pulp (ECF), white
Wrapping crepe	Pulp (ECF), white
Sewing thread	PES, white
Hand towel	2pcs, Scrim Reinforced Tissue, white

## Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

## Packaging Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PP/PE film	242 x 334
Dispenser Pack	Carton, TCF or ECF quality	2534: 381 x 279 x 320 2535: 381 x 279 x 320 2537: 381 x 279 x 320 2539: 381 x 279 x 320 2540: 381 x 279 x 370 2544: 381 x 279 x 320 2545: 381 x 279 x 320 2547: 381 x 279 x 370 2549: 381 x 279 x 370 2550: 381 x 279 x 370
Transport Carton	Carton, TCF or ECF quality	2534: 574 x 401 x 330 2535: 574 x 401 x 330 2537: 574 x 401 x 330 2539: 574 x 401 x 330 2540: 574 x 401 x 380 2544: 574 x 401 x 330 2545: 574 x 401 x 330 2547: 574 x 401 x 380 2549: 574 x 401 x 380 2550: 574 x 401 x 380
Packing Tape	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
2534	06438129112238	06438129312232	06438129512236
2535	06438129112245	06438129312249	06438129512243
2537	06438129112252	06438129312256	06438129512250
2539	06438129112269	06438129312263	06438129512267
2540	06438129112276	06438129312270	06438129512274
2544	06438129112283	06438129312287	06438129512281
2545	06438129112290	06438129312294	06438129512298
2547	06438129112306	06438129312300	06438129512304
2549	06438129112313	06438129312317	06438129512311
2550	06438129112320	06438129312324	06438129512328

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.