

evercare[®] Isolation gown Special

REF: 2733, 2734, 2735

evercare[®] Isolation gown Special is a comfortable, fluid protective gown for various uses in health care settings. The product is for single-use and non-sterile. Extra comfortable design made with knitted cuffs. Laminated in the front and sleeves for better fluid protection. Not laminated on the back to increase the breathability. Tie-bands in waist and Velcro in the neck for convenient adjustment to fit all. Available in several sizes.

Key Features

- ▶ Fluid impermeable body front and sleeves
- ▶ Cuffs around wrists
- ▶ Packed in 5 pcs per pouch
- ▶ Three sizes available
- ▶ Color: yellow



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2733	Isolation gown Special, Yellow	M, length 115 cm	5 / 25
2734	Isolation gown Special, Yellow	L, length 135 cm	5 / 25
2735	Isolation gown Special, Yellow	XL, length 145 cm	5 / 25

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



Not made with natural rubber latex



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-5
EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR Class I

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 20417, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
Body front and sleeves	Spunbond PP nonwoven, PE laminate, antistatic, yellow
Body back	Spunbond PP nonwoven, antistatic, yellow
Belt, Neck ribbon	Spunbond PP nonwoven, antistatic, yellow
Cuff	PES, white
Neck label	Coated nonwoven, white
Velcro band	PA, woven ribbon with hook or loop, white
Sewing thread	PES, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	PE bag transparent	2733: 370 x 520 x 10 2734: 330 x 460 x 10 2735: 370 x 520 x 10
Transport Carton	TCF or ECF quality	2733: 430 x 345 x 205 2734: 370 x 290 x 205 2735: 460 x 345 x 216
Packing Tape	PVC-free	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2733	N/A	06438129311327	06438129511321
2734	N/A	06438129311334	06438129511338
2735	N/A	06438129311341	06438129511345

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.