

evercare® Isolation gown Special

REF: 2733, 2734, 2735

evercare® Isolation gown Special is a comfortable, fluid protective gown for various uses in health care settings. The product is for single-use and non-sterile. Extra comfortable design made with knitted cuffs. Laminated in the front and sleeves for better fluid protection. Not laminated on the back to increase the breathability. Tie-bands in waist and Velcro in the neck for convenient adjustment to fit all. Available in several sizes.

Key Features:

- ▶ Fluid impermeable body front and sleeves
- ▶ Cuffs around wrists
- ▶ Three sizes available
- ▶ Color: yellow
- ▶ Packed in 5 pcs per pouch




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2733	Isolation gown Special, yellow	M, 115 cm	5/25
2734	Isolation gown Special, yellow	L, 135 cm	5/25
2735	Isolation gown Special, yellow	XL, 145 cm	5/25

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices



REACH regulation:

The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA)

Biological evaluation:

EN ISO 10993-5, EN ISO 10993-10

Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, Class I

Compliance with product standards:

N/A

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic
Body front & sleeves	Spunbond PP nonwoven, PE laminate, antistatic, yellow
Body back	Spunbond PP nonwoven, antistatic, yellow
Belt, Neck ribbon	Spunbond PP nonwoven, antistatic, yellow
Cuff	PES, white
Neck label	Coated nonwoven, white
Velcro band	PA, woven ribbon with hook or loop, white
Sewing thread	PES, white

Warnings and Recommendations

Warnings:

The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Storage Recommendation







Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone

Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	PE-bag	2733: 330 x 460 2734: 370 x 520 2735: 370 x 520
Transport Carton	Carton, TCF or ECF quality	2733: 370 x 290 x 205 2734: 430 x 345 x 205 2735: 460 x 345 x 216
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2733	N/A	 6 438129 311327	 6 438129 511321
2734	N/A	 6 438129 311334	 6 438129 511338
2735	N/A	 6 438129 311341	 6 438129 511345



OneMed Group Oy
 Address: Metsäläntie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
 www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1