

evercare[®] Medical face mask with tie-bands, type IIR

REF: 2884

evercare[®] Medical face mask with tie-bands, type II R provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Color: Outer facing blue, inner facing white
- ▶ Anatomic pleated folding limiting inward leakage



Tie-band



Splash resistant

| REF | DESCRIPTION | SIZE | Pack Count dispenser/carton |
|------|---------------------------------------------------------------|-------------|--------------------------------|
| 2884 | Medical face mask with tie-bands, type IIR, non-sterile, Blue | 175 x 95 mm | 50 / 1200 |

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



| | |
|---------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| REACH regulation: | The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA). |
| Biological evaluation: | EN ISO 10993-5 EN ISO 10993-10 |
| Compliance with legal requirements: | Medical Device Regulation (EU) 2017/745, MDR Class I |
| Compliance with product standards: | EN 14683:2019 |
| Quality standards: | EN ISO 13485 |
| Label information and packaging compliance: | EN ISO 15223-1, EN ISO 20417 |

Material

| Product Part | Raw Material & Characteristics |
|--------------|--------------------------------------------------------|
| Outer facing | Spunbond PP nonwoven, blue |
| Filter media | Meltblown PP nonwoven, white |
| Inner facing | Spunbond PP nonwoven, white |
| Tie-bands | Spunbond PP nonwoven, white, upper/lower length 400 mm |
| Nose wire | PP coated Galvanized iron wire, 110 mm |

Warnings & Recommendations

| | |
|--------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Warnings | The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. |
| Storage Recommendation | Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone. |
| Disposal Recommendation | Dispose in accordance with local regulations. |
| Shelf Life | 5 years |

Packing Information

| Package | Material | Size (L x W x H) mm |
|-------------------------|----------------------------|-----------------------|
| Single Pack | N/A | N/A |
| Dispenser Pack | Coated paper | 190 x 105 x 140 |
| Transport Carton | Carton, TCF or ECF quality | 780 x 335 x 300 |
| Packing Tape | PVC-free, Transparent | N/A |

Barcodes / UDI-DI

| REF | Single Pack | Dispenser | Transport Carton |
|-------------|-------------|----------------|------------------|
| 2884 | N/A | 06438129324648 | 06438129324642 |

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.