



## evercare® XP Standard OP-Gowns, sterile

REF: 2504, 2505, 2507, 2509, 2510, 2513, 2514, 2515, 2517, 2519, 2520

evercare® XP Standard OP-Gowns are suitable for dry operations. Made in a soft SMS material with antistatic and alcohol repellent treatment for enhanced user protection. All evercare® OP-Gowns provide an effective barrier to body fluids, skin particles and microbes and are designed for optimal comfort and ease of movement. The adjustable Velcro in the neck enables optimal fit while securing adjustments. High quality cuffs for improved comfort and tight fit round wrists. The gown is sterilized with EO-gas and packed in easy-to-open unit pouches. State-of-the-art packaging system with handy functional dispenser and labelling. Unit packs are supplied with barcoded removable labels for improved traceability.

### Key Features

- ▶ Good protection and Xtra comfort, antistatic and alcohol repellent treated
- ▶ Ultrasonic heat sealed seams for sleeves and shoulders to increase comfort by reducing excessive materials causing discomfort
- ▶ Gown with overlapping back and adjustable Velcro in the neck enables optimal fit for individual adjustments.
- ▶ Clear stamp with size and performance type.
- ▶ High quality cuffs for improved comfort and tight fit round wrists.
- ▶ Belt card with clear indications of sterile / non-sterile for secure donning.
- ▶ Book-folding enabling unpacking and donning without risk for contamination.
- ▶ SMS material with low linting and very high abrasion resistance




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2504	XP Standard OP-Gown, sterile, Blue	Size M-short, Length 120 cm	21 / 42
2505	XP Standard OP-Gown, sterile, Blue	Size M, Length 131 cm	20 / 40
2507	XP Standard OP-Gown, sterile, Blue	Size L, Length 142 cm	20 / 40
2509	XP Standard OP-Gown, sterile, Blue	Size XL, Length 152 cm	18 / 36
2510	XP Standard OP-Gown, sterile, Blue	Size XXL, Length 160 cm	18 / 36
2513	XP Standard OP-Gown with hand towels, sterile, Blue	Size XXL-Long, Length 170 cm	14 / 28
2514	XP Standard OP-Gown with hand towels, sterile, Blue	Size M-short, Length 120 cm	21 / 42
2515	XP Standard OP-Gown with hand towels, sterile, Blue	Size M, Length 131 cm	20 / 40
2517	XP Standard OP-Gown with hand towels, sterile, Blue	Size L, Length 142 cm	20 / 40
2519	XP Standard OP-Gown with hand towels, sterile, Blue	Size XL, Length 152 cm	18 / 36
2520	XP Standard OP-Gown with hand towels, sterile, Blue	Size XXL, Length 160 cm	18 / 36

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class Is
Compliance with product standards:	EN ISO 11135, EN 13795-1 Standard performance
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

## Material

Product Part	Raw Material & Characteristics
Body, sleeve, belt	PP SMS nonwoven, blue
Cuff	PES, rib knitted fabric
Neck ribbon	PP Spunbond nonwoven, white
Velcro band	PA woven ribbon with hook or loop
Belt card	Cellulose pulp (ECF), white
Wrapping crepe	Pulp (ECF), white
Sewing thread	PES, white
Hand towel	2pcs, Scrim Reinforced Tissue, white

## Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

## Packaging Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PP/PE film	242 x 334
Dispenser Pack	Carton, TCF or ECF quality	381 x 279 x 320
Transport Carton	Carton, TCF or ECF quality	574 x 401 x 330
Packing Tape	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
2504	06438129112139	06438129312133	06438129512137
2505	06438129112146	06438129312140	06438129512144
2507	06438129112153	06438129312157	06438129512151
2509	06438129112160	06438129312164	06438129512168
2510	06438129112177	06438129312171	06438129512175
2513	06438129124293	06438129324297	06438129524291
2514	06438129112184	06438129312188	06438129512182
2515	06438129112191	06438129312195	06438129512199
2517	06438129112207	06438129312201	06438129512205
2519	06438129112214	06438129312218	06438129512212
2520	06438129112221	06438129312225	06438129512229

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.