

439

Non-woven swabs

Compressas em não tecido

Compresses en non tissé

Vliesstoff-Kompressen

Nonwoven kompres

Compressa in T.N.T.

Compresa de tela sin tejer

Non-woven kompress



Mais imagens



Non-woven swabs

- For wound dressing and cleansing. For treatment of lightly or moderately exuding wounds.
- 70% viscose, 30% polyester.
- Approx. weight: 30, 40 and 50g/sq.m.
- Some models with X-Ray detectable thread with approx. 60% barium sulfate contents.
- Some models available with Y shape cut.
- Paper bag + transport carton.

Banded by 5

Swabs banded by 5 are available for ETO sterilization only.
Specific code numbers and quotations will apply.





Code	Size	Plies	g/m2	X-Ray	Y shape cut	Pcs./Bag	Pcs./Tr. carton
439-001	5 x 5cm	4	30			100	5000
439-002	7,5 x 7,5cm	4	30			100	5000
439-003	7,5 x 10cm	4	30			100	5000
439-004	10 x 10cm	4	30			100	5000
439-005	10 x 20cm	4	30			100	5000
439-014	10 x 10cm	8	30			100	1000
439-050	5 x 5cm	4	40			100	5000
439-051	7,5 x 7,5cm	4	40			100	5000
439-052	7,5 x 10cm	4	40			100	5000
439-053	10 x 10cm	4	40			100	5000
439-054	10 x 20cm	4	40			100	5000
439-060	5 x 5cm	8	40			100	1000
439-061	7,5 x 7,5cm	8	40			100	1000
439-063	10 x 10cm	8	40			100	1000
439-074	10 x 20cm	6	40			100	1000
439-100	5 x 5cm	4	50			100	5000
439-101	7,5 x 7,5cm	4	50			100	5000
439-102	7,5 x 10cm	4	50			100	5000
439-103	10 x 10cm	4	50			100	5000
439-104	10 x 20cm	4	50			100	5000
439-151	7,5 x 7,5cm	4	50	X		100	5000
439-153	10 x 10cm	4	50	X		100	5000
439-154	10 x 20cm	4	50	X		100	5000
439-163	10 x 10cm	8	50	X		100	1000
439-173	15 x 15cm	4	30			100	3000
439-174	15 x 20cm	4	30			100	2500
439-175	15 x 25cm	4	30			100	3200
439-176	20 x 20cm	4	30			100	2000
439-177	20 x 25cm	4	30			100	1800
439-178	25 x 25cm	4	30			100	1500
439-180	15 x 15cm	6	30			100	2000
439-181	15 x 20cm	6	30			100	2000
439-182	15 x 25cm	6	30			100	2400
439-183	20 x 20cm	6	30			100	1800
439-184	20 x 25cm	6	30			100	1800
439-188	15 x 20cm	4	30	X		100	2500
439-190	20 x 25cm	4	30	X		100	1800
439-202	7,5 x 7,5cm	4	30		X	100	5000
439-204	10 x 10cm	4	30		X	100	5000
439-213	20 x 20 cm	4	40			100	2000