

evercare[®] Examination Gloves, Vinyl SOFT

REF: 6300, 6301, 6302, 6303, 6304

Soft powder-free vinyl examination glove for single use.

Key Features

- ▶ Soft vinyl glove
- ▶ Ambidextrous and soft vinyl glove with rolled rim for extra cuff strength
- ▶ Powder free and latex free
- ▶ Good alternative to those suffering from type I and accelerator related type IV allergies



Recommendations for use








Basic cleaning



Basic procedures
and wards



Wash

REF	DESCRIPTION	SIZE CODE	Overall Length	Pack Count dispenser/carton
6300	Examination Gloves, Vinyl SOFT, XS, Clear		240 mm	100 / 1000
6301	Examination Gloves, Vinyl SOFT, S, Clear		240 mm	100 / 1000
6302	Examination Gloves, Vinyl SOFT, M, Clear		240 mm	100 / 1000
6303	Examination Gloves, Vinyl SOFT, L, Clear		240 mm	100 / 1000
6304	Examination Gloves, Vinyl SOFT, XL, Clear		240 mm	100 / 1000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Freedom of holes:	AQL 1.5
Force at break (N):	≥ 3.6 N (before and after challenge test)
Biological evaluation:	Total quantity of powder residues for powder-free glove: < 2 mg/glove Cytotoxicity tested acc. to EN ISO 10993-5 Irritation and skin sensitization tested acc. to EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 455: 1 – 4
Food safety:	The glove is suitable for handling aqueous foodstuffs - not for handling fatty foods.
Other tests performed:	Viral penetration according to ASTM F 1671
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Palm	Polyvinyl Chloride, single wall thickness min. 0.08 mm
Finger	Polyvinyl Chloride, single wall thickness min. 0.08 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Cardboard, TCF or ECF quality	240 x 120 x 65
Transport Carton	Corrugated board, TCF or ECF quality	250 x 250 x 340
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
6300	N/A	06438129310757	06438129510751
6301	N/A	06438129310764	06438129510768
6302	N/A	06438129310771	06438129510775
6303	N/A	06438129310788	06438129510782
6304	N/A	06438129310795	06438129510799

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.