

evercare® Examination Gloves, Vinyl SENSITIVE X-LONG

REF: 6305, 6306, 6307, 6308, 6309

Vinyl x-long powder-free and phthalate-free examination glove for single use.

Key Features

- ▶ Long
- ▶ Ambidextrous with rolled rim for extra cuff strength
- ▶ Phthalate-free, powder-free
- ▶ Latex-free
- ▶ Good alternative to those suffering from type I or accelerator relating type IV allergy



Recommendations for use








Basic cleaning



Basic procedures
and wards



Wash

REF	DESCRIPTION	SIZE CODE	Overall Length	Pack Count dispenser/carton
6305	Examination Gloves, Vinyl SENSITIVE X-LONG, XS, Clear		290 mm	100 / 1000
6306	Examination Gloves, Vinyl SENSITIVE X-LONG, S, Clear		290 mm	100 / 1000
6307	Examination Gloves, Vinyl SENSITIVE X-LONG, M, Clear		290 mm	100 / 1000
6308	Examination Gloves, Vinyl SENSITIVE X-LONG, L, Clear		290 mm	100 / 1000
6309	Examination Gloves, Vinyl SENSITIVE X-LONG, XL, Clear		290 mm	100 / 1000

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.







EN 455

REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Freedom of holes:	AQL 1.5
Force at break (N):	≥ 3.6 N (before and after challenge test)
Biological evaluation:	Total quantity of powder residues for powder-free glove: < 2 mg/glove Cytotoxicity tested acc. to EN ISO 10993-5 Irritation and skin sensitization tested acc. to EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 455: 1 – 4
Food safety:	In accordance with EU Regulations 1935/2004 and 10/2011 the glove is suitable for handling all foodstuffs except fatty foods.
Other tests performed:	Viral penetration according to ASTM F 1671
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Palm	Polyvinyl Chloride, single wall thickness min. 0.07 mm
Finger	Polyvinyl Chloride, single wall thickness min. 0.08 mm

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Cardboard, TCF or ECF quality	260 x 110 x 85
Transport Carton	Corrugated board, TCF or ECF quality	439 x 230 x 270
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
6305	N/A	06438129310801	06438129510805
6306	N/A	06438129310818	06438129510812
6307	N/A	06438129310825	06438129510829
6308	N/A	06438129310832	06438129510836
6309	N/A	06438129310849	06438129510843

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.