

# evercare® Examination Gloves, Nitrile SVEA

REF: 6415, 6416, 6417, 6418, 6419

Nitrile examination glove and protective glove for single use, blue/lilac.

## Key Features

- ▶ Powder-free, ambidextrous and with rolled rim for extra cuff strength.
- ▶ Excellent tactile sensitivity and enhanced wet and dry grip with micro-textured finger tips.
- ▶ Excellent fit, easy to don and remove.
- ▶ Well recognized examination glove color for ward.
- ▶ Higher quality with AQL 1.0 of pinholes for safer use.
- ▶ Very good viral barrier properties, protects against a wide variety of chemicals.
- ▶ Free of latex and proteins, good alternative to those suffering from type I allergies.



## Recommendations for use



Critical care



Dental care



Medical laboratory








Pharmacy



Cleaning and surface disinfection



Good protection against infection

REF	DESCRIPTION	SIZE CODE	OVERALL LENGTH	PACK COUNT dispenser/carton
6415	Examination Gloves, Nitrile SVEA, XS, Blue/Lilac		240 mm	150 / 1500
6416	Examination Gloves, Nitrile SVEA, S, Blue/Lilac		240 mm	150 / 1500
6417	Examination Gloves, Nitrile SVEA, M, Blue/Lilac		240 mm	150 / 1500
6418	Examination Gloves, Nitrile SVEA, L, Blue/Lilac		240 mm	150 / 1500
6419	Examination Gloves, Nitrile SVEA, XL, Blue/Lilac		240 mm	150 / 1500

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices and personal protective equipments.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Freedom of holes:</b>	AQL 1.0
<b>Force at break (N):</b>	≥ 6.0 N (before and after challenge test)
<b>Biological evaluation:</b>	Total quantity of powder residues for powder-free glove: < 2 mg/glove Cytotoxicity tested acc. to EN ISO 10993-5 Irritation and skin sensitization tested acc. to EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I Personal Protective Equipment Regulation (EU) 2016/425, PPE R category III EU Declaration of Conformity for Personal Protective Equipment: <a href="https://www.evercare-medical.com/about-us">https://www.evercare-medical.com/about-us</a>
<b>Compliance with product standards:</b>	EN 455: 1 – 4
<b>Chemical resistance:</b>	Tested in accordance with EN 16523-1
<b>Food safety:</b>	In accordance with EU Regulation 1935/2004 the glove is suitable for handling all foodstuffs.
<b>Other tests performed:</b>	Air leak and Water leak according to EN 374-2 (Level 2) Chemicals residue Dexterity, length and pH according to EN 420 Mechanical performance according to EN 388 Microbiological quality RoHS compliance Resistance to degradation by chemicals according to EN 374-4 Viral penetration according to EN 374-5 and ASTM F1671
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041, EN 420

## Material

Product Part	Raw Material & Characteristics
<b>Palm</b>	Acrylonitrile-butadiene (Nitrile), single wall thickness min. 0.07 mm
<b>Finger</b>	Acrylonitrile-butadiene (Nitrile), single wall thickness min. 0.10 mm

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. The product contains dithiocarbamates – not to be used in case of hypersensitivity to these substances.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A
Dispenser Pack	Cardboard, TCF or ECF quality	240 x 120 x 75
Transport Carton	Corrugated board, TCF or ECF quality	385 x 250 x 250
Packing Tape	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
6415	N/A	6438129313109	6438129513103
6416	N/A	6438129313116	6438129513110
6417	N/A	6438129313123	6438129513127
6418	N/A	6438129313130	6438129513134
6419	N/A	6438129313147	6438129513141

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.