

429

Gauze balls

Bolas de gazeTampons de gazeMulltupferGaasdepperTampone tondoBola de gasaGasv vstork- boll

Gauze balls

- Cotton, elastic ring.
- For fluid absorption during surgery, cleaning and disinfection of the skin.
- 17 threads/sq.cm with approx. 23g/sq.m,
- Raw edges of the fabric are turned inside and kept with the help of an elastic ring made of silicone.
- Plastic bag + transport carton.

DIAMETER:
DI  METRO:



Flat size.
Di  metro achatado.



Round size.
Di  metro arredondado.

Non sterile

Code	REF	Weight	X Diamater	Y Diamenter	Threads	Pcs./Bag	Pcs./Tr. carton
429-004	D	0,37	25	35	17	100	1000
429-005	E	0,71	30	40	17	100	1000
429-006	F	0,83	40	50	17	100	1000
429-007	G	1,63	50	60	17	100	1000
429-008	H	1,87	55	70	17	100	1000
429-009	I	2,92	65	80	17	100	1000
429-010	J	3,32	75	90	17	100	1000
429-011	K	3,78	85	100	17	100	1000