

SELEFA® Gauze swabs, sterile

REF: 2222702, 2222705, 2222722, 2222725, 2222742, 2222745, 2222762, 2222765, 222294, 222295, 222296, 222298, 222299, 222300

SELEFA® Gauze swabs, sterile are primarily intended to cover or clean and absorb blood and exudate from wounds and abrasions. Can also be used for non-invasive medical procedures, such as wound padding, drying, protection, skin and wound cleansing or as secondary dressing. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze.

Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ 17 threads per cm²
- ▶ Folded edges to prevent loose threads
- ▶ Sterilized using EO
- ▶ Easy-to-use unit pouch




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2222702	Gauze swabs 2-pack, 8-ply, sterile, White	5x5 cm	120 / 2400
2222705	Gauze swabs 5-pack, 8-ply, sterile, White	5x5 cm	125 / 2500
2222722	Gauze swabs 2-pack, 8-ply, sterile, White	7,5x7,5 cm	120 / 2400
2222725	Gauze swabs 5-pack, 8-ply, sterile, White	7,5x7,5 cm	125 / 2500
2222742	Gauze swabs 2-pack, 8-ply, sterile, White	10x10 cm	120 / 2160
2222745	Gauze swabs 5-pack, 8-ply, sterile, White	10x10 cm	125 / 2250
2222762	Gauze swabs 2-pack, 12-ply, sterile, White	10x20 cm	40 / 720
2222765	Gauze swabs 5-pack, 12-ply, sterile, White	10x20 cm	40 / 720
222294	Gauze swabs, 5-pack, 12-ply, sterile, White	10x10 cm	80 / 1280
222295	Gauze swabs, 5-pack, 12-ply, sterile, White	7,5X7,5CM	80 / 1600
222296	Gauze swabs, 5-pack, 12-ply, sterile, White	5x5 cm	80 / 1600
222298	Gauze swabs, 5-pack, 8-ply, sterile, White	20x80 cm	25 / 150
222299	Gauze swabs, 5-pack, 8-ply, sterile, White	10x60 cm	65 / 390
222300	Gauze swabs, 25-pack, 8-ply, sterile, White	10x10 cm	150 / 2400

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR Class Is
Compliance with product standards:	EN 14079, EN ISO 11135, EN ISO 11737-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Main material	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	2222702: 75 x 90
		2222705: 75 x 90
		2222722: 100 x 114
		2222725: 100 x 114
		2222742: 133 x 135
		2222745: 133 x 135
		2222762: 133 x 241
		2222765: 133 x 241
		222294: 133 x 135
		222295: 100 x 114
		222296: 75 x 90
		222298: 133 x 240
		222299: 133 x 240
		222300: 133 x 135
Dispenser Pack	Carton, TCF or ECF quality	2222702: 90 x 78 x 145
		2222705: 90 x 78 x 145
		2222722: 115 x 95 x 145
		2222725: 115 x 95 x 145
		2222742: 135 x 130 x 145
		2222745: 135 x 130 x 145
		2222762: 250 x 130 x 90
		2222765: 250 x 130 x 90
		222294: 135 x 130 x 145
		222295: 115 x 95 x 145
		222296: 90 x 78 x 145
		222298: 240 x 130 x 210
		222299: 240 x 130 x 210
		222300: 140 x 135 x 160
Transport Carton	Carton, TCF or ECF quality	2222702: 400 x 190 x 300
		2222705: 400 x 190 x 300
		2222722: 600 x 200 x 300
		2222725: 600 x 200 x 300
		2222742: 400 x 300 x 420
		2222745: 400 x 300 x 420
		2222762: 400 x 280 x 510
		2222765: 400 x 280 x 510
		222294: 540 x 320 x 500
		222295: 540 x 320 x 500
		222296: 410 x 250 x 450
		222298: 400 x 250 x 430
		222299: 400 x 250 x 430
		222300: 550 x 330 x 290
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
2222702	06438129150971	06438129350975	06438129550979
2222705	06438129151015	06438129351019	06438129551013
2222722	06438129150988	06438129350982	06438129550986
2222725	06438129151022	06438129351026	06438129551020
2222742	06438129150995	06438129350999	06438129550993
2222745	06438129151039	06438129351033	06438129551037
2222762	06438129151008	06438129351002	06438129551006
2222765	06438129151046	06438129351040	06438129551044
222294	06438129151442	06438129351446	06438129551440
222295	06438129151459	06438129351453	06438129551457
222296	06438129151466	06438129351460	06438129551464
222298	06438129151473	06438129351477	06438129551471
222299	06438129151480	06438129351484	06438129551488
222300	06438129151497	06438129351491	06438129551495

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.