



## SELEFA® Nonwoven swabs, 2-pack, sterile

REF: 2232702, 2232722, 2232742, 2232762

SELEFA® Nonwoven swabs, 4-ply, sterile are produced from soft and highly absorptive nonwoven material. They are primarily used for cleaning and exudate absorption of wounds and abrasions.

### Key Features

- ▶ High quality viscose/polyester nonwoven material
- ▶ Low linting and high softness
- ▶ Keeps its shape, even when wet
- ▶ Free of binding agents and optical brighteners
- ▶ Effective in liquid and exudate absorption
- ▶ Less adherent to the wound than gauze devices
- ▶ Available in several sizes
- ▶ For wound cleansing and padding, exudate absorption, pre-operative skin cleansing etc.
- ▶ Can be used as a primary or secondary dressings
- ▶ Sterilized using EO
- ▶ Easy-to-use unit pouch
- ▶ 2-pack



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2232702	Nonwoven swabs, 4-ply, 2-pack, sterile, White	5x5 cm	140 / 2800
2232722	Nonwoven swabs, 4-ply, 2-pack, sterile, White	7,5x7,5 cm	140 / 2800
2232742	Nonwoven swabs, 4-ply, 2-pack, sterile, White	10x10 cm	140 / 2520
2232762	Nonwoven swabs, 4-ply, 2-pack, sterile, White	10x20 cm	70 / 1260

Sales Unit: Transport Carton

## Regulations

N/A	
<div> <div> <div>CE</div> <div>2460</div> </div> <div> <div>CE</div> </div> <div> <div>STERILE</div> <div>EO</div> </div> <div> <div>MD</div> </div> <div> <div>LATEX</div> <div>Not made with natural rubber latex</div> </div> <div> <div></div> </div> <div> <div></div> </div> <div> <div></div> </div> <div> <div>15 °C</div> <div>25 °C</div> </div> <div> <div>2</div> <div>STERILE</div> </div> <div> <div></div> </div> </div>	
REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN ISO 11607-1, EN ISO 11135
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, N/A

## Material

Product Part	Raw Material & Characteristics
Main material	Nonwoven, 70% viscose and 30% polyester, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	2232702: 98 x 116 2232722: 72 x 88 2232742: 98 x 113 2232762: 122 x 138
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	2232702: 237 x 100 x 116 2232722: 90 x 70 x 180 2232742: 115 x 95 x 180 2232762: 140 x 120 x 180
<b>Transport Carton</b>	Carton, TCF or ECF quality	2232702: 595 x 245 x 415 2232722: 370 x 190 x 370 2232742: 595 x 200 x 370 2232762: 435 x 370 x 370
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>2232702</b>	06438129151053	06438129351057	06438129551051
<b>2232722</b>	06438129151060	06438129351064	06438129551068
<b>2232742</b>	N/A	N/A	N/A
<b>2232762</b>	06438129151084	06438129351088	06438129551082

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.