



SELEFA® Surgical Tape PE-film

REF: 2402031, 2402032, 2402033

SELEFA® Surgical Tape PE-film is an economical, general purpose surgical tape. It can be used for fixation of tubes, catheters and small medical equipment as well as for fixation of all kinds of dressings. It is breathable and can easily be removed without any adhesive rests.

Key Features

- ▶ Latex-free
- ▶ Breathable
- ▶ Excellent adhesive properties
- ▶ Thin transparent perforated PE-film




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
2402031	Surgical Tape PE-film, Transparent	1,25 cm x 9,1 m	24 / 720
2402032	Surgical Tape PE-film, Transparent	2,5 cm x 9,1 m	12 / 360
2402033	Surgical Tape PE-film, Transparent	5 cm x 9,1 m	6 / 180

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	EN 1939, EN ISO 14184-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Tape	LLDPE film, semi-transparent
Glue	Acrylic pressure sensitive adhesive, latex free
Core	PE, cream white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	158 x 105 x 54
Transport Carton	Carton, TCF or ECF quality	335 x 335 x 285
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2402031	N/A	06438129352115	06438129552119
2402032	N/A	06438129352122	06438129552126
2402033	N/A	06438129352139	06438129552133

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.