

SELEFA® Surgical Tape Nonwoven

REF: 2402021, 2402022, 2402023, 2405021, 2405022, 2405023

SELEFA® Surgical tape Nonwoven is an economical, general purpose surgical tape. It can be used for fixation of tubes, catheters and small medical equipment as well as for fixation of all kinds of dressings. It is breathable and can easily be removed without any adhesive rests.

Key Features

- ▶ Latex-free
- ▶ Breathable
- ▶ Excellent adhesive properties
- ▶ Thin nonwoven material
- ▶ Available with or without holder




| REF | DESCRIPTION | SIZE | Pack Count dispenser/carton |
|---------|---|-----------------|--------------------------------|
| 2402021 | Surgical Tape Nonwoven, White | 1,25 cm x 9,1 m | 24 / 720 |
| 2402022 | Surgical Tape Nonwoven, White | 2,5 cm x 9,1 m | 12 / 360 |
| 2402023 | Surgical Tape Nonwoven, White | 5 cm x 9,1 m | 6 / 180 |
| 2405021 | Surgical Tape Nonwoven with holder, White | 1,25 cm x 9,1 m | 24 / 720 |
| 2405022 | Surgical Tape Nonwoven with holder, White | 2,5 cm x 9,1 m | 12 / 360 |
| 2405023 | Surgical Tape Nonwoven with holder, White | 5 cm x 9,1 m | 6 / 180 |

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



| | |
|---|---|
| REACH regulation: | The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA). |
| Biological evaluation: | EN ISO 10993-5 EN ISO 10993-10 |
| Compliance with legal requirements: | Medical Device Regulation (EU) 2017/745, MDR class I |
| Compliance with product standards: | EN 1939, EN ISO 14184-1 |
| Quality standards: | EN ISO 13485 |
| Label information and packaging compliance: | EN ISO 15223-1, EN ISO 20417 |

Material

| Product Part | Raw Material & Characteristics |
|--------------|---|
| Tape | Paper and wet-laid nonwoven, white |
| Glue | Acrylic pressure sensitive adhesive, latex free |
| Core | PE, cream white |
| Holder | PP, blue |

Warnings & Recommendations

| | |
|--------------------------------|---|
| Warnings | The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. |
| Storage Recommendation | Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone. |
| Disposal Recommendation | Dispose in accordance with local regulations. |
| Shelf Life | 3 years |

Packing Information

| Package | Material | Size (L x W x H) mm |
|-------------------------|----------------------------|--|
| Single Pack | N/A | N/A |
| Dispenser Pack | Carton, TCF or ECF quality | 2402021: 136 x 90 x 54 2402022: 136 x 90 x 54 2402023: 136 x 90 x 54 2405021: 145 x 142 x 78 2405022: 145 x 142 x 66 2405023: 145 x 142 x 60 |
| Transport Carton | Carton, TCF or ECF quality | 2402021: 290 x 290 x 285 2402022: 290 x 290 x 285 2402023: 290 x 290 x 285 2405021: 445 x 310 x 410 2405022: 445 x 310 x 350 2405023: 445 x 310 x 320 |
| Packing Tape | PVC-free, Transparent | N/A |

Barcodes / UDI-DI

| REF | Single Pack | Dispenser | Transport Carton |
|---------|-------------|----------------|------------------|
| 2402021 | N/A | 06438129352085 | 06438129552089 |
| 2402022 | N/A | 06438129352092 | 06438129552096 |
| 2402023 | N/A | 06438129352108 | 06438129552102 |
| 2405021 | N/A | 06438129352146 | 06438129552140 |
| 2405022 | N/A | 06438129352153 | 06438129552157 |
| 2405023 | N/A | 06438129352160 | 06438129552164 |

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.