

## evercare<sup>®</sup> Elastic fixation bandage Fix Premium

REF: 4110, 4111, 4112, 4113, 4114, 4150, 4151, 4152, 4153, 4154

evercare<sup>®</sup> Elastic fixation bandage Fix Premium is stretchable with soft and porous structure. Due to its air permeability and softness the bandages are easy to apply, conformable and comfortable to use. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. Package designed for easy product identification.

### Key Features

- ▶ Fully woven edges to minimize loose threads
- ▶ Fixation of all types of dressings, especially for joints and round body parts
- ▶ Soft and comfortable
- ▶ Non-slip fit
- ▶ 100% stretchability
- ▶ 60% polyamide and 40% viscose




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4110	Elastic fixation bandage Fix Premium, polyamide and viscose, single pack, non-sterile, White	4 cm x 4 m	20 / 720
4111	Elastic fixation bandage Fix Premium, polyamide and viscose, single pack, non-sterile, White	6 cm x 4 m	20 / 560
4112	Elastic fixation bandage Fix Premium, polyamide and viscose, single pack, non-sterile, White	8 cm x 4 m	20 / 480
4113	Elastic fixation bandage Fix Premium, polyamide and viscose, single pack, non-sterile, White	10 cm x 4 m	20 / 400
4114	Elastic fixation bandage Fix Premium, polyamide and viscose, single pack, non-sterile, White	12 cm x 4 m	20 / 320
4150	Elastic fixation bandage Fix Premium, polyamide and viscose, White	4 cm x 4 m	20 / 720
4151	Elastic fixation bandage Fix Premium, polyamide and viscose, White	6 cm x 4 m	20 / 560
4152	Elastic fixation bandage Fix Premium, polyamide and viscose, White	8 cm x 4 m	20 / 480
4153	Elastic fixation bandage Fix Premium, polyamide and viscose, White	10 cm x 4 m	20 / 400
4154	Elastic fixation bandage Fix Premium, polyamide and viscose, White	12 cm x 4 m	20 / 320

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-5 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	DIN 61632, DIN 61634
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
<b>Bandage</b>	60 % polyamide, 40 % viscose, Unstretched length 2 m, Stretched length 4 m, Elongation approx. 100%, Elasticity approx. 85%

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4110: 175 x 140 x 45 4111: 175 x 140 x 65 4112: 175 x 140 x 85 4113: 175 x 140 x 105 4114: 175 x 140 x 125 4150: 175 x 140 x 45 4151: 175 x 140 x 65 4152: 175 x 140 x 85 4153: 175 x 140 x 105 4154: 175 x 140 x 125
<b>Transport Carton</b>	Carton, TCF or ECF quality	4110: 425 x 365 x 295 4111: 475 x 365 x 295 4112: 530 x 365 x 295 4113: 545 x 365 x 295 4114: 520 x 365 x 295 4150: 425 x 365 x 295 4151: 475 x 365 x 295 4152: 530 x 365 x 295 4153: 545 x 365 x 295 4154: 520 x 365 x 295
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4110	N/A	06438129314564	06438129514568
4111	N/A	06438129314571	06438129514575
4112	N/A	06438129314588	06438129514582
4113	N/A	06438129314595	06438129514599
4114	N/A	06438129314601	06438129514605
4150	N/A	06438129314632	06438129514636
4151	N/A	06438129314649	06438129514643
4152	N/A	06438129314656	06438129514650
4153	N/A	06438129314663	06438129514667
4154	N/A	06438129314670	06438129514674

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.