

## evercare<sup>®</sup> Elastic fixation bandage Selexif, sterile

REF: 4130-01, 4131-01, 4132-01, 4133-01, 4134-01, 4135-01

evercare<sup>®</sup> Elastic fixation bandage Selexif, sterile is stretchable with soft and porous structure. Due to its air permeability and softness the bandages are easy to apply, conformable and comfortable to use. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. Package designed for easy product identification.

### Key Features

- ▶ Fully woven edges to minimize loose threads
- ▶ Fixation of all types of dressings, especially for joints and round body parts
- ▶ Soft and comfortable
- ▶ Non-slip fit
- ▶ 90% stretchability
- ▶ 100% polyester
- ▶ Sterilized with EO
- ▶ Single, peel open packed




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4130-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	4 cm x 4 m	25 / 100
4131-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	6 cm x 4 m	25 / 100
4132-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	8 cm x 4 m	25 / 100
4133-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	10 cm x 4 m	25 / 100
4134-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	12 cm x 4 m	25 / 100
4135-01	Elastic fixation bandage Selexif, polyester, single pack, sterile, White	15 cm x 4 m	10 / 60

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5  
EN ISO 10993-10

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class Is

**Compliance with product standards:** DIN 61632, DIN 61634, EN 1049-2

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
<b>Bandage</b>	100% polyester, white, Unstretched length 2 m, Stretched length 4 m, Elongation approx. 90%, Elasticity approx. 80%

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	4130-01: 57 x 95
		4131-01: 57 x 118
		4132-01: 57 x 143
		4133-01: 57 x 143
		4134-01: 57 x 188
		4135-01: 57 x 188
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4130-01: 190 x 90 x 160
		4131-01: 190 x 110 x 160
		4132-01: 190 x 140 x 160
		4133-01: 190 x 140 x 160
		4134-01: 190 x 190 x 160
		4135-01: 190 x 190 x 160
<b>Transport Carton</b>	Carton, TCF or ECF quality	4130-01: 335 x 205 x 195
		4131-01: 335 x 205 x 235
		4132-01: 335 x 205 x 295
		4133-01: 335 x 205 x 295
		4134-01: 395 x 205 x 335
		4135-01: 395 x 205 x 335
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4130-01	06438129121896	06438129321890	06438129521894
4131-01	06438129121902	06438129321906	06438129521900
4132-01	06438129121919	06438129321913	06438129521917
4133-01	06438129121926	06438129321920	06438129521924
4134-01	06438129121933	06438129321937	06438129521931
4135-01	6438129121940	6438129321944	6438129521948

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.