

## evercare<sup>®</sup> Elastic fixation bandage Selexif

REF: 4181-01, 4182-01, 4183-01, 4184-01, 4185-01

evercare<sup>®</sup> Elastic fixation bandage Selexif is stretchable with soft and porous structure. Due to its air permeability and softness the bandages are easy to apply, conformable and comfortable to use. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. Package designed for easy product identification.

### Key Features

- ▶ Fully woven edges to minimize loose threads
- ▶ Fixation of all types of dressings, especially for joints and round body parts
- ▶ Soft and comfortable
- ▶ Non-slip fit
- ▶ 90% stretchability
- ▶ 100% polyester



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4181-01	Elastic fixation bandage Selexif, polyester, non-sterile, White	4 cm x 4 m	20 / 600
4182-01	Elastic fixation bandage Selexif, polyester, non-sterile, White	6 cm x 4 m	20 / 560
4183-01	Elastic fixation bandage Selexif, polyester, non-sterile, White	8 cm x 4 m	20 / 480
4184-01	Elastic fixation bandage Selexif, polyester, non-sterile, White	10 cm x 4 m	20 / 400
4185-01	Elastic fixation bandage Selexif, polyester, non-sterile, White	12 cm x 4 m	20 / 320

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-5 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	DIN 61632, DIN 61634, EN 1049-2
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
<b>Bandage</b>	100% polyester, white, Unstretched length 2 m, Stretched length 4 m, Elongation approx. 100%, Elasticity approx. 80%

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4181-01: 165 x 45 x 135 4182-01: 165 x 65 x 135 4183-01: 165 x 85 x 135 4184-01: 165 x 105 x 135 4185-01: 165 x 125 x 135
<b>Transport Carton</b>	Carton, TCF or ECF quality	4181-01: 425 x 350 x 245 4182-01: 485 x 350 x 290 4183-01: 540 x 350 x 290 4184-01: 550 x 350 x 290 4185-01: 520 x 350 x 290
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>4181-01</b>	N/A	06438129321333	06438129521337
<b>4182-01</b>	N/A	06438129321340	06438129521344
<b>4183-01</b>	N/A	06438129321357	06438129521351
<b>4184-01</b>	N/A	06438129321364	06438129521368
<b>4185-01</b>	N/A	06438129321371	06438129521375

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.