

evercare[®] Self-adhesive fabric

REF: 4555, 4560, 4565, 4570, 4575

evercare[®] Self-adhesive fabric is used to support fixation of dressings, swabs, catheters, tubes and cannulae. Skin-friendly adhesive, stretchable and conforms well to body contours. Allows movement when in place. Easy to remove release paper premeasured in 1 cm intervals. Breathable, pleasant to wear and moisture repellent.

Key Features

- ▶ Ideal for areas such as neck, elbow and knee
- ▶ Perforated split for easy removal and application
- ▶ Barrier for blood and fluid strike through



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4555	Self-adhesive fabric roll, White	2,5 cm x 10 m	1 / 120
4560	Self-adhesive fabric roll, White	5 cm x 10 m	1 / 60
4565	Self-adhesive fabric roll, White	10 cm x 10 m	1 / 36
4570	Self-adhesive fabric roll, White	15 cm x 10 m	1 / 24
4575	Self-adhesive fabric roll, White	20 cm x 10 m	1 / 18

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



Not made with natural rubber latex



REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-5
EN ISO 10993-10

Compliance with legal requirements:

Medical Device Regulation (EU) 2017/745, MDR class I

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Carrier	100% polyester, nonwoven spunlance, water repellent
Coating	Pressure sensitive glue, thermosetting acrylic solution polymer, latex free
Release paper	Coated paper, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	4555: 30 x 83 x 83 4560: 55 x 83 x 83 4565: 105 x 83 x 83 4570: 155 x 83 x 83 4575: 205 x 83 x 83
Transport Carton	Carton, TCF or ECF quality	4555: 360 x 330 x 280 4560: 360 x 300 x 280 4565: 360 x 330 x 280 4570: 360 x 330 x 280 4575: 435 x 275 x 280
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4555	N/A	06438129313833	06438129513837
4560	N/A	06438129313840	06438129513844
4565	N/A	06438129313857	06438129513851
4570	N/A	06438129313864	06438129513868
4575	N/A	06438129313871	06438129513875

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.