


evercare® Extremity set

REF: 1231-02

evercare® Extremity set, sterile is a sterile draping set for knee operations.

Key Features

- ▶ Elastic fenestration
- ▶ Absorbent patch
- ▶ Elastic stockinette
- ▶ Tube holders
- ▶ Adhesive OP Tapes
- ▶ Main drape complies to EN 13795-1 High performance critical product area requirements



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1231-02	Extremity set, sterile	See Material table	4 / 8

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-01 EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN13795-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
4 Hand towels	30 x 34 cm, Pulp, natural white, thread reinforced
1 Mayo cover	78 x 145 cm, reinforced, Laminate: blue PP nonwoven/blue-white PE film/blue PE film, blue absorbent patch 65 x 85 cm
2 Adhesive Op-tapes	9 x 50 cm, Nonwoven polyester, blue, adhesive coated
1 Instrument table cover, folded	150 x 190 cm, PP nonwoven/polyethylene film laminate, blue
1 Stockinette small	22 x 75 cm, Styrene Butadiene Styrene film/elastic gauze polyester, teleskop folding
1 Extremity drape	230 x 300 cm, PP nonwoven/polyethylene film laminate, blue/PP nonwoven, white, blue absorbent patch 50 x 100 cm, PP nonwoven, elastic aperture 7 cm, elastic material styrene/butadiene/styrene, tube holders in blue nonwoven
1 Instrument table cover	150 x 190 cm, PP nonwoven/polyethylene film laminate, blue

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	420 x 603
Dispenser Pack	TCF or ECF quality	485 x 285 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1231-02	06438129116779	06438129316773	06438129516777

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.