



## SELEFA® Optimia Gauze squares

REF: 222301, 222302, 222303, 222304, 222305

SELEFA® Optimia Gauze squares are primarily intended to cover or clean and absorb blood and exudate from wounds and abrasions. Can also be used for non-invasive medical procedures, such as wound padding, drying, protection, skin and wound cleansing or as secondary dressing. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze.

### Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ Unfolded
- ▶ Can be sterlized using EO
- ▶ Excelent absorption
- ▶ 20 threads per cm<sup>2</sup>




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
222301	Optimia Gauze squares, White	10x30 cm	500 / 10000
222302	Optimia Gauze squares, White	20x30 cm	500 / 10000
222303	Optimia Gauze squares, White	30x30 cm	500 / 5000
222304	Optimia Gauze squares, White	40x30 cm	500 / 4000
222305	Optimia Gauze squares, White	50x50 cm	250 / 2500

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-5 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	EN 14079
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
Main material	100 % cotton gauze, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Medical grade paper, White	222301: 150 x 150 x 100 222302: 200 x 200 x 100 222303: 300 x 300 x 100 222304: 400 x 400 x 100 222305: 500 x 500 x 100
<b>Transport Carton</b>	Carton, TCF or ECF quality	222301: 340 x 340 x 500 222302: 430 x 430 x 500 222303: 630 x 320 x 500 222304: 830 x 420 x 420 222305: 520 x 520 x 500
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>222301</b>	N/A	06438129354188	06438129554182
<b>222302</b>	N/A	06438129354195	06438129554199
<b>222303</b>	N/A	06438129354201	06438129554205
<b>222304</b>	N/A	06438129354218	06438129554212
<b>222305</b>	N/A	06438129354225	06438129554229

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.