

evercare[®] Transparent IV film dressings, reinforced, sterile

REF: 4740

evercare[®] Comfort I.V. dressing, reinforced, sterile is designed for wound treatment or cannula fixation. Adhesive film dressing maintains optimal moisture level of treated area while protecting against dust, water and bacteria at the same time. Polyurethane(PU) semi-permeable foil used for its composition allows the patient to take a shower reducing the risk of wound contamination. The dressing's transparency allows visual inspection of the wound without removal. The dressing design ensures painless removal. The transparent foil material is reinforced with nonwoven material at one side to strengthen the fixation.

Key Features

- ▶ Soft and flexible
- ▶ High breathability
- ▶ Acrylic based adhesive
- ▶ Easy-to-use unit pouch
- ▶ Effective barrier to exudate, water and bacteria
- ▶ Primary dressing for cannula fixation
- ▶ Contain extra strips for enhanced securement over the I.V. site
- ▶ Reinforced with split
- ▶ High MVTR factor




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4740	Transparent IV film dressings, reinforced, sterile, Transparent	7x9 cm	100 / 1200

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10
EN ISO 10993-11

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class IIa

Compliance with product standards: EN ISO 11135, EN ISO 11737-1, EN 13726-2, EN 13726-4

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Carrier	Transparent PU film, breathable
Coating	Solvent based acrylic adhesive, transparent, breathable
Reinforced carrier part	PET nonwoven, white
Coating reinforcement	Solvent based acrylic adhesive, transparent, breathable
Release paper	Medical grade paper silicone coated, white
Backing liner release foil	PE/PET or PP film, transparent
Date label	Non-woven fabric with adhesive, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	155 x 98
Dispenser Pack	Carton, TCF or ECF quality	105 x 100 x 160
Transport Carton	Carton, TCF or ECF quality	350 x 250 x 330
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4740	06438129124354	06438129324358	06438129524352

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.