



SELEFA® Apron without sleeves, FELIX

REF: 210435

SELEFA® Apron without sleeves, FELIX is a disposable apron made of fluid resistant polyethylene that help to protect clothing from spills and splashes during medical procedures and cleaning tasks. Also ideal for use in food preparations and in industrial applications.

Key Features

- ▶ Disposable fluid resistant Apron
- ▶ Convenient pullover
- ▶ Transparent, Embossed
- ▶ Available on a roll of 50




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
210435	Apron without sleeves, FELIX, embossed, Transparent	81 x 150 cm	50 / 500

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for personal protective equipments.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Compliance with legal requirements:	Personal Protective Equipment Regulation (EU) 2016/425, PPE R category I. EU Declaration of Conformity for Personal Protective Equipment: https://www.evercare-medical.com/about-us
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
Main material	LLDPE with LDPE, weight 18.60 g/m ²

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations. Clean and dry product can be recycled among PE-LD plastics.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Corrugated cardboard	430 x 90 x 90
Transport Carton	Corrugated cardboard	476 x 435 x 194
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
210435	N/A	6438129350180	6438129550184

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.